


Armrest Assembly



Pos.	Leckey Code	Description	Qty	Compatible with	
				Size 1	Size 2
1	117-715	Armrest Moulded Assembly	1	✓	✓
2	137-9011	Armrest Adjustment Kit	1	✓	✓
3	137-9013	Armrest Locking Kit	1	✓	✓
4	137-9012	Armrest Outer Adjustment Assembly	1	✓	✓
5	SJ109	1/4 Single Popper	1	✓	✓
6	117-145	Amrest Outer Extrusion - Size 1	1	✓	
6	137-146	Amrest Outer Extrusion - Size 2	1		✓
7	117-9005	Amrrest Adjuster - Size 1	1	✓	
7	137-9014	Amrrest Adjuster - Size 2	1		✓