


Hi-Low Chassis - Powered



Pos.	Leckey Code	Description	Qty	Compatible with	
				Size 1	Size 2
1	102-1693	Chassis Welded Assy - Size 1 & 2	1	✓	✓
2	IA13	Castors	4	✓	✓
3	102-1694	Actuator Assy	1	✓	✓
4	SI74	Charger - UK	1	✓	✓
5	SI125	Charger - EUR	1	✓	✓
6	SI126	Charger - USA/Japan	1	✓	✓
7	SI127	Charger - AUS	1	✓	✓
8	102-1695	500N Blocolift - Size 1 & 2 Assy	1	✓	✓
9	102-1696	800N Blocolift Size 3 Assy	1	✓	✓
10	SI64	Battery	1	✓	✓
11	SI67	Handset	1	✓	✓
12	117-651	Lever Arm Assembly Size 1-3	1	✓	✓
13	117-633	Mygo Interface Assembly	1	✓	✓