


Dynamic Backrest (50N/100N/150N)



Pos.	Leckey Code	Description	Qty	Compatible with	
				Size 1	Size 2
1	137-660	Dynamic Backrest Assy (50N Gas Spring)	1	✓	✓
1	137-661	Dynamic Backrest Assy (100N Gas Spring)	1	✓	✓
1	137-662	Dynamic Backrest Assy (150N Gas Spring)	1	✓	✓
2	SP-1060	Dynamic Backrest Assy Attachment Kit	1	✓	✓
3	SP-1061	Dynamic Backrest Adjustment Kit	1	✓	✓
4	SP-1004	Back Angle Attachment Kit	1	✓	✓