


## Chassis Assembly



Pos.	Leckey Code	Description	Qty	Compatible with	
				Size 1	Size 2
1	IA13	BLACK "LECKEY" TANG	4	✓	✓
2	SF56	1-1/2" RIBBED TUBE	4	✓	✓
3	SH62	1 1/4" ROUND BLACK	4	✓	✓
4	143-9000	CHASSIS ASSEMBLY - SIZE 1	1	✓	
4	143-9002	CHASSIS ASSEMBLY - SIZE 2	1		✓
5	143-9001	CHASSIS FIXINGS KIT	1	✓	✓