


Comfee Backrests



Pos.	Leckey Code	Description	Qty	Compatible with		
				Size 1	Size 2	Size 3
1	149-1949	S1 Backrest One Piece Kit	1	✓		
2	149-1950	S1 Backrest Two Piece Kit	1	✓	✓	
3	149-2949	S2 Backrest One Piece Kit	1		✓	
4	149-2950	S2 Backrest Two Piece Kit	1			✓
4	149-3949	S3 Backrest One Piece Kit	1			✓
4	149-3950	S3 Backrest Two Piece Kit	1			✓
5	149-1951	Handle Bar Assembly	1	✓	✓	✓
5	149-1952	Headrest Clamp Kit	1	✓	✓	✓
5	149-1953	Backrest Height Clamp Kit	1	✓	✓	✓
6	149-2800	Gas Spring Cable for Backrest Angle	1	✓	✓	✓
6	SI161	Mechlok Cable for Tilt in Space	1	✓	✓	✓