

LECKEY[®]

Easy Seat



LECKEY^o

Stronger together.

We work with individuals, therapists and carers to design products with both a clinical and an emotional focus. Using the latest research and clinical understanding, we create practical solutions which are easily integrated into family life, because for us, life is about going, enjoying, participating and doing.



24 hour
postural care for
babies, kids & adults.
Sleeping, Sitting,
Standing, Walking,
Moving, Bathing,
Toileting.



The Easy Seat is a robust, stable and adjustable seat for children with mild to moderate needs, aged 1-11 years of age.

Available in 4 sizes, the Easy Seat is available with either fixed or castor legs.

The product is designed with seat depth & seat width adjustment, back height & angle adjustment and accommodation for lower leg length growth.

A large range of accessories is available including flexible & rigid laterals, side pads, pommel, headrest, footplate and sandal straps ensuring that the product meets a wide range of clinical needs.





Backrest with height and angle adjustment

Wipe clean cushions

Height adjustable armrest

Side pads

Seat depth adjustment

Robust steel frame

45° pelvic belt

Height adjustable legs

Footplate and sandals



Angle adjustable backrest



Additional headrest option.



Footplate and sandals



45° pelvic belt



Height adjustable legs



Height adjustable castor legs



Lightweight activity tray

The Easy Seat was specifically designed to provide an attractive seating solution without compromising the postural needs of the less involved child.

The high level of adjustability and large range of accessories available for the Easy Seat allow it to grow with the child.

The robust and sturdy Easy Seat design is available in four sizes from nursery to school age children. An important feature is that it should grow with the child.

Therefore, as standard, we have incorporated seat depth adjustment, lower leg adjustment and backrest height adjustment.

The height adjustability of the Easy Seat not only accommodates lower leg lengths but allows integration within the nursery environment and, as the footplate remains close to the floor, eases self transfer.

The Easy Seat is available in four sizes.

The standard Easy Seat includes:
Adjustable seat and backrest cushions;
45° pelvic belt; side pads; fixed or castor legs.
Sizes 1 & 2 come as standard in Meadow fabric.
Sizes 3 & 4 come in Mondeo Blue.

Product Sizing

Size	1		2		3		4	
Age	1 - 3		2 - 5		3 - 8		6 - 11	
Codes	LES/1/FL LES/1/CL		LES/2/FL LES/2/CL		LES/3/FL LES/3/CL		LES/4/FL LES/4/CL	
	Fixed Legs	Castor Legs						
	mm	inches	mm	inches	mm	inches	mm	inches
Easy Seat								
Length	400	16	480	19	570	22	620	24
Width	350	14	370	14	450	18	490	19
Weight of Product	12kg		15kg		20kg		22kg	
User Standing Height								
min	760	30	930	37	1000	40	1150	46
max	1000	40	1150	46	1300	52	1450	58
Max User Weight								
	30kg	66lbs	30kg	66lbs	40kg	88lbs	40kg	88bs
Backrest Height (sitting shoulder height)								
min	225	9	265	10	325	13	400	16
max	285	11	340	14	425	17	500	20
Seat Base Depth (backrest to back of knee)								
min	175	7	225	9	250	10	300	12
max	250	10	300	12	350	14	400	16
Seat Base Width (hip width)								
min	125	5	150	6	225	9	275	11
max	225	9	250	10	325	13	375	15
Calf Length (sole of foot to back of knee)								
min	150	6	220	9	250	10	300	12
max	230	9	300	12	350	14	420	17
Tray Height (seat to top of tray)								
min	135	5	135	5	180	7	180	7
max	225	9	225	9	325	13	325	13
Armrest Height (seat base to elbow)								
max	110	4	120	5	160	6	180	7
min	150	6	180	7	210	8	240	10
Degree of angular adjustment on backrest								
Prone	5°		5°		5°		5°	
Recline	10°		10°		10°		10°	
Castor size								
	75	3	75	3	75	3	75	3
Tray sizes								
	450 x420	18 x16	450 x420	18 x16	590 x520	23 x20	590 x520	23 x20

Leckey Servicing

Our FREE service and support includes:
Free product training
Free product set up
Free product assessment
Free product re-assessment
Free repair within 2 year warranty

For service policies on all products outside warranty, please contact Leckey's Customer Service department. James Leckey Design Ltd as manufacturer with sole responsibility declares that all products conform to 93/42/EEC guidelines and EN12182 technical aids for disabled persons general requirements and test methods. Order forms and spare parts lists to extend the service life of the product and allow reissue are available on request or online at www.leckey.com.

Accessories



LES/#/A Height adjustable tray



LES/#/C Adjustable pommel



LES/#/D Flexible lateral supports



LES/#/E Rigid lateral supports



LES/#/F Footplate
LES/#/N With sandal straps



LES/#/I Grab rail



LES/#/L Flat high back
LES/#/M Laterals for Flat high back



LES/#/P Grab posts



LES/#/Q Push handle



LES/#/S Potty insert

Fabric Options



Meadow



Mondeo Blue



LECKEY^o


Established in 1983, Leckey is a globally recognised pioneer in the research and development of products that help adults and children with disabilities to go, do, enjoy and participate in everyday activities throughout the day and night.

We take a highly clinical approach to product design and development. Through in-depth clinical research studies with leading universities, and extensive trials with occupational therapists, physiotherapists, users and their families, we continue to develop posturally supportive, family friendly products for all day care, at every stage of life.

Through early intervention, childhood and adulthood Leckey's experienced team of designers, therapists and bioengineers work together to develop products that meet the clinical needs of the healthcare professionals and the social needs of the user.

To achieve this, we work with the healthcare professionals, the individuals and carers who use our products everyday. With their help, we create the dependable, durable, proven and high performance products that we are known for worldwide.





24 hour postural
care for babies,
kids & adults.
Sleeping, Sitting,
Standing, Walking,
Moving, Bathing,
Toileting.