


## Foot Support Assembly



Pos.	Leckey Code	Description	Qty	Compatible with	
				Size 1	Size 2
1	133-2071	Sandal Assy	1	✓	✓
2	133-2072	Sandal Bolt Kit	1	✓	✓
3	133-2073	Footplate Adjustment Clamp Kit	1	✓	✓
4	133-2074	Inner Leg LH	1	✓	✓
4	133-2075	Inner Leg RH	1	✓	✓
5	133-2076	Leg Height Adjustment Clamp Kit	1	✓	✓
6	133-654	Outer Leg	1	✓	✓