


PAL -Tray Assembly



Pos.	Leckey Code	Description	Qty	Compatible with		
				Size 1	Size 2	Size 3
1	135-1620	Activity Tray Size 1	1	✓		
1	135-2620	Activity Tray Size 2	1		✓	
1	135-3620	Activity Tray Size 3	1			✓
2	SF90	Wixroyd Knob M6	1	✓	✓	✓
3	135-1656	Tray Pad Size 1	1	✓		
3	135-2656	Tray Pad Size 2	1		✓	
3	135-3656	Tray Pad Size 3	1			✓