

Whitmyer Sizing Guide

Whitmyer Adjust-A-Plush

Adjust-A-PLUSH Standard: P8088L Narrow: P8044L	<p>Standard: 50 – 152 mm Narrow: 44 – 100 mm</p>	Standard Adjustment: 200 mm – 460 mm 127 mm	Narrow Adjustment: 152 – 355 mm 76 mm
---	--	--	--

Whitmyer Plush

480 mm PLUSH P1900L	<p>95 mm</p>	<p>178 mm 250 mm</p>	<p>228 mm 190 mm 482 mm</p>	<p>44 mm</p>
350 mm PLUSH P1400L	<p>63 mm</p>	<p>100 mm 190 mm</p>	<p>127 mm 336 mm</p>	<p>44 mm</p>
250 mm PLUSH P1000L	<p>32 mm</p>	<p>152 mm</p>	<p>127 mm 235 mm</p>	<p>44 mm</p>

Whitmyer Contoured Cradle

Large P400L	<p>108 mm</p>	<p>100 mm 171 mm</p>	<p>176 mm 190 mm 330 mm</p>	<p>25 mm</p>
Medium P300L	<p>82 mm</p>	<p>82 mm 152 mm</p>	<p>120 mm 158 mm 250 mm</p>	<p>15 mm</p>
Small P200L	<p>57 mm</p>	<p>82 mm 127 mm</p>	<p>100 mm 139 mm 216 mm</p>	<p>15 mm</p>

Whitmyer Heads Up

Heads Up Pads				
Large				
Short				
Long / Tail				
Medium				
Short				
Long / Tail				
Small				
Short				
Long / Tail				
Early Intervention				
Short				
Long / Tail				



Product Compatibility Guide

<u>Product</u>	<u>Compatibility</u>
Seating	
JCM Triton	Compatible using 8 Inch PRO Bar
JCM Jupiter	Not Compatible
JCM Neptune	Compatible using 8 Inch PRO Bar
JCM Sunbeam	Not Compatible
JCM Star	Not Compatible
Leckey Mygo +	Compatible using 8 Inch PRO Bar
Leckey Mygo	Compatible using 8 Inch PRO Bar
Leckey Mygo Max	Compatible using 8 Inch PRO Bar
Leckey Everyday Activity Seat	Compatible using 8 Inch PRO Bar
Leckey Squiggles Seat	Compatible using 4/8 Inch PRO Bar
Standing	
Leckey Squiggles +	Compatible using 4 Inch PRO Bar
Leckey Squiggles TT	Compatible using 4 Inch PRO Bar
Leckey Mygo Stander	Compatible using 4/8 Inch PRO Bar
Leckey Horizon Stander	Compatible using 8 Inch PRO Bar
JCM Aurora	Compatible using AXYS Mounting System
Leckey Squiggles Stander	Compatible using 4 Inch PRO Bar