



Chassis Footprints

User Instructions



CHASSIS FOOTPRINTS DIMENSIONS

Product	Size	A (mm)	B (mm)	A (Inches)	B (Inches)
Squiggles Seat on Hi Low Chassis	N/A	540	640	21.3	25.2
Mygo Seat on High Low Chassis	All	570	710	22.4	28
Mygo + on High Low Chassis	All	570	710	22.4	28
Everyday Seat	1	570	780	22.4	30.7
Everyday Seat	2	570	780	22.4	30.7
Everyday Seat	3	660	890	26	35
KIT Seat on Hydraulic Chassis	1 and 2	610	690	24	27.2
KIT Seat on Powered Chassis	1 and 2	610	730	24	28.7
Mygo Max Seat on Standard Chassis	N/A	835	1245	32.9	49
Mygo Max Seat on KIT Hydraulic Chassis	N/A	610	690	24	27.2
Mygo Max Seat on KIT Powered Chassis	N/A	610	730	24	28.7
Comfee Seat [Discontinued]	1	710	870	28	34.3
Comfee Seat [Discontinued]	2	760	945	29.9	37.2
Comfee Seat [Discontinued]	3	785	1000	30.9	39.4
Freestander	1	650	650	25.6	25.6
Freestander	2	750	750	29.5	29.5
Freestander	3	750	750	29.5	29.5
Pronestander	1	60	820	2.4	32.3
Pronestander	2	650	1070	25.6	42.1
Pronestander	3	700	1270	27.6	50
Squiggles Stander	N/A	550	800	21.7	31.5
Squiggles + Stander	N/A	580	960	22.8	37.8
Squiggles TT	N/A	580	960	22.8	37.8
Mygo Stander	1	630	978	24.8	38.5
Mygo Stander	2	650	1130	25.6	44.5
Horizon Stander	1	580	1135	22.8	44.7
Horizon Stander	2	640	1280	25.2	50.4
Horizon Stander	3	640	1280	25.2	50.4
MyWay Walker	1	570	720	22.4	28.3
MyWay Walker	2	670	850	26.4	33.5
MyWay Walker	3	720	1100	28.3	43.3





Leckey,
19c Ballinderry Road
Lisburn, BT28 2SA
Northern Ireland

leckey.com

{+44} 28 9260 0750
hello@leckey.com

