



Horizon Derotational Pelvic Band

User Instructions



HORIZON DEROTATIONAL PELVIC BAND SIZES

Code	DESCRIPTION	Dimensions	Min Hip Circumference	Max Hip Circumference
LHZ/1/181-XXX	SUPINE PELVIC BAND SIZE 1			
	HIP STRAP OUTER	1530 x 100		
	HIP STRAP INNER	990 x 100	500 [20"]	800 [32"]
LHZ/1/182-XXX	PRONE PELVIC BAND 1			
	PRONE HIP STRAP OUTER FABRIC	1530 x 100		
	PRONE HIP STRAP INNER FABRIC	990 x 100	500 [20"]	800 [32"]
LHZ/2/181-XXX	SUPINE PELVIC BAND SIZE 2			
	HIP STRAP OUTER	1680 x 100		
	HIP STRAP INNER	1230 x 100	750 [30"]	950 [38"]
LHZ/2/182-XXX	PRONE PELVIC BAND 2			
	PRONE HIP STRAP OUTER FABRIC	1680 x 100		
	PRONE HIP STRAP INNER FABRIC	1230 x 100	750 [30"]	950 [38"]
LHZ/3/181-XXX	SUPINE PELVIC BAND SIZE 3			
	HIP STRAP OUTER	1680 x 100		
	HIP STRAP INNER	1230 x 100	1000 [40"]	950 [38"]
LHZ/3/182-XXX	PRONE PELVIC BAND 3			
	PRONE HIP STRAP OUTER FABRIC	1680 x 100		
	PRONE HIP STRAP INNER FABRIC	1230 x 100	1000 [40"]	950 [38"]

*XXX replaced by colour

GB (Green)
BB (Blue)
PB (Pink)
GRB (Grey)





Leckey,
19c Ballinderry Road
Lisburn, BT28 2SA
Northern Ireland

leckey.com

{+44} 28 9260 0750
hello@leckey.com

