



Leckey Advanced Bath Chair

Name	
Date:	
Email	
Contact Number	
Invoice Address	
Delivery Address	

Quote Order

Purchase Order	
----------------	--



	Item Code	Description	Quantity	Price
<input type="checkbox"/>	LSF/1	Leckey Shower Frame (for Advance Bath Chair only)		XXX
Frame				
<input type="checkbox"/>	LAB/1/FR	Advanced Bath Chair Frame Size 1		XXX
<input type="checkbox"/>	LAB/2/FR	Advanced Bath Chair Frame Size 2		XXX
<input type="checkbox"/>	LAB/3/FR	Advanced Bath Chair Frame Size 3		XXX
<input type="checkbox"/>	LAB/4/FR	Advanced Bath Chair Frame Size 4		XXX
Cover				
<input type="checkbox"/>	LAB/1/D-PK	Advanced Bath Chair Size 1 Cover - Pink		XXX
<input type="checkbox"/>	LAB/2/D-RB	Advanced Bath Chair Size 2 Cover - Royal Blue		XXX
<input type="checkbox"/>	LAB/2/D-PK	Advanced Bath Chair Size 2 Cover - Pink		XXX
<input type="checkbox"/>	LAB/3/D-RB	Advanced Bath Chair Size 3 Cover - Royal Blue		XXX
<input type="checkbox"/>	LAB/3/D-PK	Advanced Bath Chair Size 3 Cover - Pink		XXX
<input type="checkbox"/>	LAB/4/D-RB	Advanced Bath Chair Size 4 Cover - Royal Blue		XXX
<input type="checkbox"/>	LAB/4/D-PK	Advanced Bath Chair Size 4 Cover - Pink		XXX
Hip Belt				
<input type="checkbox"/>	LAB/1/E-RB	Hip Belt - Royal Blue		XXX
<input type="checkbox"/>	LAB/1/E-PK	Hip Belt - Pink		XXX
Chest Belt				
<input type="checkbox"/>	LAB/1/A-RB	Advanced Bath Chair Chest Belt - Royal Blue		XXX
<input type="checkbox"/>	LAB/1/A-PK	Advanced Bath Chair Chest Belt - Pink		XXX
Head Lateral Support				
<input type="checkbox"/>	LAB/1/C-RB	Advanced Bath Chair Head Laterals for Advance Bath Chair Royal Blue		XXX
<input type="checkbox"/>	LAB/1/C-PK	Advanced Bath Chair Head Laterals for Advance Bath Chair Pink		XXX
Pommel/Abduction Guides				
<input type="checkbox"/>	LAB/1/B-RB	Advanced Bath Chair Size 1 & 2 Abduction Belt for Advance Bath Chair Royal Blue		XXX
<input type="checkbox"/>	LAB/1/B-PK	Advanced Bath Chair Size 1 & 2 Abduction Belt for Advance Bath Chair Pink		XXX
<input type="checkbox"/>	LAB/2/B-RB	Advanced Bath Chair Size 3 & 4 Abduction Belt for Advance Bath Chair Royal Blue		XXX
<input type="checkbox"/>	LAB/2/B-PK	Advanced Bath Chair Size 3 & 4 Abduction Belt for Advance Bath Chair Pink		XXX
Package				
<input type="checkbox"/>	LAB/1/COMP-RB	Size 1 Advance Bath Chair Including Hip and Chest Strap - Royal Blue		XXX
<input type="checkbox"/>	LAB/1/COMP-PK	Size 1 Advance Bath Chair Including Hip and Chest Strap - Pink		XXX
<input type="checkbox"/>	LAB/2/COMP-RB	Size 2 Advance Bath Chair Including Hip and Chest Strap - Royal Blue		XXX
<input type="checkbox"/>	LAB/2/COMP-PK	Size 2 Advance Bath Chair Including Hip and Chest Strap - Pink		XXX
<input type="checkbox"/>	LAB/3/COMP-RB	Size 3 Advance Bath Chair Including Hip and Chest Strap - Royal Blue		XXX
<input type="checkbox"/>	LAB/3/COMP-PK	Size 3 Advance Bath Chair Including Hip and Chest Strap - Pink		XXX
<input type="checkbox"/>	LAB/4/COMP-RB	Size 4 Advance Bath Chair Including Hip and Chest Strap - Royal Blue		XXX
<input type="checkbox"/>	LAB/4/COMP-PK	Size 4 Advance Bath Chair Including Hip and Chest Strap - Pink		XXX

Technical Specifications

Specification	Size 1	Size 2	Size 3	Size 4
Age	1 - 5 years	4 - 9 years	8 - 14 years	12 - 18 years
User Height (Max / Min)	750 mm / 1050 mm	950 mm / 1250 mm	1150 mm / 1450 mm	1350 mm / 1650 mm
Head Support	150 mm	150 mm	215 mm	215 mm
Back Support Height	350 mm	400 mm	460 mm	550 mm
Seat Depth	225 mm	275 mm	335 mm	400 mm
Calf Support	260 mm	260 mm	340 mm	340 mm
Seat Height	50 mm - 425 mm	50 mm - 425 mm	50 mm - 425 mm	50 mm - 425 mm
Maximum Weight	72 kgs	72 kgs	72 kgs	72 kgs
Total Width	430 mm	430 mm	430 mm	430 mm