



Leckey Mygo Stander

Name	
Date:	
Email	
Contact Number	
Invoice Address	
Delivery Address	

Quote Order



	Item Code	Description	Quantity	Price
Frame				
<input type="checkbox"/>	143-1600	Mygo Stander size 1 stander - with manual angle adjustment PU chest pad PU hip pad PU kr		XXX
<input type="checkbox"/>	143-2600	Mygo Stander size 2 stander - with manual angle adjustment PU chest pad PU hip pad PU k		XXX
Tray				
<input type="checkbox"/>	143-1633	Mygo Stander tray with bowl insert - size 1		XXX
<input type="checkbox"/>	143-1609	Mygo Stander Mygo stander tray - size 1		XXX
<input type="checkbox"/>	143-2633	Mygo Stander Mygo stander tray with bowl insert - size 2		XXX
<input type="checkbox"/>	143-2609	Mygo Stander Mygo stander tray - size 2		XXX
<input type="checkbox"/>	143-1649	Mygo Stander grab rail		XXX
Covers				
<input type="checkbox"/>	143-1770-08	Mygo Stander Size 1 Cover pack - blue		XXX
<input type="checkbox"/>	143-1770-07	Mygo Stander Size 1 Cover pack - orange		XXX
<input type="checkbox"/>	143-1770-09	Mygo Stander Size 1 Cover pack - pink		XXX
<input type="checkbox"/>	143-2770-08	Mygo Stander Size 2 Cover pack - blue		XXX
<input type="checkbox"/>	143-2770-07	Mygo Stander Size 2 Cover pack - orange		XXX
<input type="checkbox"/>	143-2770-09	Mygo Stander Size 2 Cover pack - pink		XXX
Chest Support				
<input type="checkbox"/>	143-1619	Mygo Stander size 1 chest lateral hardware only (pair)		XXX
<input type="checkbox"/>	143-2619	Mygo Stander size 2 chest lateral hardware only (pair)		XXX
<input type="checkbox"/>	143-1675	Mygo Stander Flipaway laterals hardware (pair)		XXX
<input type="checkbox"/>	143-1674	Mygo Stander Complex lateral vertical - left (single inc cover)		XXX
<input type="checkbox"/>	143-1673	Mygo Stander Complex lateral vertical - right (single inc cover)		XXX
<input type="checkbox"/>	143-1780-08	Mygo Stander size 1 chest/hip wraparound harness - blue		XXX
<input type="checkbox"/>	143-1780-07	Mygo Stander size 1 chest/hip wraparound harness - orange		XXX
<input type="checkbox"/>	143-1775-14	Mygo Stander size 1 chest/hip wraparound harness - Black Infection Control		XXX
<input type="checkbox"/>	143-1780-09	Mygo Stander size 1 chest/hip wraparound harness - pink		XXX
<input type="checkbox"/>	143-2780-08	Mygo Stander size 2 chest/hip wraparound harness - blue		XXX
<input type="checkbox"/>	143-2780-07	Mygo Stander size 2 chest/hip wraparound harness - orange		XXX
<input type="checkbox"/>	143-2775-14	Mygo Stander size 2 chest/hip wraparound harness - Black Infection Control		XXX
<input type="checkbox"/>	143-2780-09	Mygo Stander size 2 chest/hip wraparound harness - pink		XXX
<input type="checkbox"/>	143-3780-08	Mygo Stander size 3 chest/hip wraparound harness - Blue		XXX
<input type="checkbox"/>	143-3780-07	Mygo Stander size 3 chest/hip wraparound harness - Orange		XXX
<input type="checkbox"/>	143-3775-14	Mygo Stander size 3 chest/hip wraparound harness - Black Infection Control		XXX
<input type="checkbox"/>	143-3780-09	Mygo Stander size 3 chest/hip wraparound harness - Pink		XXX
Hip Lateral				
<input type="checkbox"/>	143-1663	Mygo Stander hip lateral hardware only		XXX
<input type="checkbox"/>	143-1780-08	Mygo Stander size 1 chest/hip wraparound harness - blue		XXX
<input type="checkbox"/>	143-1780-07	Mygo Stander size 1 chest/hip wraparound harness - orange		XXX
<input type="checkbox"/>	143-1775-14	Mygo Stander size 1 chest/hip wraparound harness - Black Infection Control		XXX
<input type="checkbox"/>	143-1780-09	Mygo Stander size 1 chest/hip wraparound harness - pink		XXX
<input type="checkbox"/>	143-2780-08	Mygo Stander size 2 chest/hip wraparound harness - blue		XXX
<input type="checkbox"/>	143-2780-07	Mygo Stander size 2 chest/hip wraparound harness - orange		XXX
<input type="checkbox"/>	143-2775-14	Mygo Stander size 2 chest/hip wraparound harness - Black Infection Control		XXX
<input type="checkbox"/>	143-2780-09	Mygo Stander size 2 chest/hip wraparound harness - pink		XXX
<input type="checkbox"/>	143-3780-08	Mygo Stander size 3 chest/hip wraparound harness - Blue		XXX
<input type="checkbox"/>	143-3780-07	Mygo Stander size 3 chest/hip wraparound harness - Orange		XXX
<input type="checkbox"/>	143-3775-14	Mygo Stander size 3 chest/hip wraparound harness - Black Infection Control		XXX
<input type="checkbox"/>	143-3780-09	Mygo Stander size 3 chest/hip wraparound harness - Pink		XXX
Derotational Belts				
<input type="checkbox"/>	143-1798-08	Mygo Stander size 1 prone derotation hip Band - blue		XXX
<input type="checkbox"/>	143-1798-07	Mygo Stander size 1 prone derotation hip Band - orange		XXX
<input type="checkbox"/>	143-1798-09	Mygo Stander size 1 prone derotation hip Band - pink		XXX
<input type="checkbox"/>	143-2798-08	Mygo Stander size 2 prone derotation hip Band - blue		XXX

<input type="checkbox"/>	143-2798-07	Mygo Stander size 2 prone derotation hip Band - orange		XXX
<input type="checkbox"/>	143-2798-09	Mygo Stander size 2 prone derotation hip Band - pink		XXX
<input type="checkbox"/>	143-3798-08	Mygo Stander size 3 prone derotation hip Band - blue		XXX
<input type="checkbox"/>	143-3798-07	Mygo Stander size 3 prone derotation hip Band - orange		XXX
<input type="checkbox"/>	143-3798-09	Mygo Stander size 3 prone derotation hip Band - pink		XXX
<input type="checkbox"/>	143-1783-08	Mygo Stander size 1 supine derotation hip Band - blue		XXX
<input type="checkbox"/>	143-1783-07	Mygo Stander size 1 supine derotation hip Band - orange		XXX
<input type="checkbox"/>	143-1783-09	Mygo Stander size 1 supine derotation hip Band - pink		XXX
<input type="checkbox"/>	143-2783-08	Mygo Stander size 2 supine derotation hip Band - blue		XXX
<input type="checkbox"/>	143-2783-07	Mygo Stander size 2 supine derotation hip Band - orange		XXX
<input type="checkbox"/>	143-2783-09	Mygo Stander size 2 supine derotation hip Band - pink		XXX
<input type="checkbox"/>	143-3783-08	Mygo Stander size 3 supine derotation hip Band - blue		XXX
<input type="checkbox"/>	143-3783-07	Mygo Stander size 3 supine derotation hip Band - orange		XXX
<input type="checkbox"/>	143-3783-09	Mygo Stander size 3 supine derotation hip Band - pink		XXX
Knee Pads				
<input type="checkbox"/>	143-1776-14	Mygo Stander size 1 basic knee straps - infection control (black)		XXX
<input type="checkbox"/>	143-1776-08	Mygo Stander size 1 basic knee straps - blue		XXX
<input type="checkbox"/>	143-1776-07	Mygo Stander size 1 basic knee straps - orange		XXX
<input type="checkbox"/>	143-1776-09	Mygo Stander size 1 basic knee straps - pink		XXX
<input type="checkbox"/>	143-2776-14	Mygo Stander size 2 basic knee straps - infection control (black)		XXX
<input type="checkbox"/>	143-2776-08	Mygo Stander size 2 basic knee straps - blue		XXX
<input type="checkbox"/>	143-2776-07	Mygo Stander size 2 basic knee straps - orange		XXX
<input type="checkbox"/>	143-2776-09	Mygo Stander size 2 basic knee straps - pink		XXX
<input type="checkbox"/>	143-1784-08	Mygo Stander split knee straps - blue - fits both sizes		XXX
<input type="checkbox"/>	143-1784-07	Mygo Stander split knee straps - orange - fits both sizes		XXX
<input type="checkbox"/>	143-1784-09	Mygo Stander split knee straps - pink - fits both sizes		XXX
<input type="checkbox"/>	143-1690	Mygo Stander Extended Knee brackets		XXX
Hip Infill				
<input type="checkbox"/>	143-1666	Mygo Stander Hip Infill Pad hardware only - requires cover		XXX
<input type="checkbox"/>	143-1799-08	Mygo Stander Size 1 Blue Hip Infill Pad cover only		XXX
<input type="checkbox"/>	143-1799-07	Mygo Stander Size 1 Orange Hip Infill Pad cover only		XXX
<input type="checkbox"/>	143-1799-09	Mygo Stander Size 1 Pink Hip Infill Pad cover only		XXX
<input type="checkbox"/>	143-2799-08	Mygo Stander Size 2 Blue Hip Infill Pad cover only		XXX
<input type="checkbox"/>	143-2799-07	Mygo Stander Size 2 Orange Hip Infill Pad cover only		XXX
<input type="checkbox"/>	143-2799-09	Mygo Stander Size 2 Pink Hip Infill Pad cover only		XXX
Shoulder Pad				
<input type="checkbox"/>	143-1638	Mygo Stander shoulder pad		XXX
<input type="checkbox"/>	143-1768-08	Mygo Stander shoulder pad cover - blue		XXX
<input type="checkbox"/>	143-1768-07	Mygo Stander shoulder pad cover - orange		XXX
<input type="checkbox"/>	143-1768-09	Mygo Stander shoulder pad cover - pink		XXX
Harness				
<input type="checkbox"/>	155-N001-L600	Small Butterfly Harness & attachment straps		XXX
<input type="checkbox"/>	155-N001-L601	Medium Butterfly Harness & attachment straps		XXX
<input type="checkbox"/>	155-N001-L602	Large Butterfly Harness & attachment straps		XXX
<input type="checkbox"/>	155-N001-L603	X-Large Butterfly Harness & attachment straps		XXX
Foot Support				
<input type="checkbox"/>	152-1600	Sandals Size 1		XXX
<input type="checkbox"/>	152-2600	Sandals Size 2		XXX
<input type="checkbox"/>	152-3600	Sandals Size 3		XXX
<input type="checkbox"/>	LHZ/1/100	Size 1 open-back sandals		XXX
<input type="checkbox"/>	LHZ/2/100	Size 2 open-back sandals		XXX
<input type="checkbox"/>	LHZ/3/100	Size 3 open-back sandals		XXX
<input type="checkbox"/>	AH01	Ankle Huggers		XXX
<input type="checkbox"/>	173-2305	Small Ankle Wraps (use with sandals)		XXX
<input type="checkbox"/>	173-3305	Large Ankle Wraps (use with sandals)		XXX
<input type="checkbox"/>	143-1615	Sandal Raiser for Mygo Stander		XXX
Headrest				
<input type="checkbox"/>	137-699-05	Contoured headrest Assembly and Cushion		XXX
<input type="checkbox"/>	137-699-05VIN	Contoured headrest Assembly and Cushion - Vinyl		XXX
<input type="checkbox"/>	117-699-05	Contoured Occipital Headrest		XXX
<input type="checkbox"/>	117-699-05VIN	Contoured Occipital Headrest - Vinyl		XXX
<input type="checkbox"/>	117-851-05	Flat Headrest Assembly and Cushion		XXX
<input type="checkbox"/>	117-851-05VIN	Flat Headrest Assembly and Cushion - Vinyl		XXX
<input type="checkbox"/>	120-893-05	Flat Headrest Laterals with covers		XXX
<input type="checkbox"/>	120-893-05VIN	Flat Headrest Laterals with covers Vinyl		XXX
<input type="checkbox"/>	117-449	1" Wooden spacer for 8 Inch PRO Bar		XXX

	Size 1	Size 2
Age	4-10	8-14
User Height	1080mm-1450mm (42.5"-57")	1350mm-1700mm (53.1"-66.9")
Max. User Weight	50kg (110lbs)	60kg (132lbs)
Chest Support Height (footplate to armpit)	760mm-1080mm (29.9"-42.5")	980mm-1300mm (38.5"-51.1")
Chest Width (distance between laterals)	160mm-280mm (6.2"-11")	200mm-320mm (7.8"-12.5")
Hip Support Height (footplate to mid-buttocks)	460mm-730mm (18.1"-28.7")	620mm-875mm (24.4"-34.4")
Hip Width (distance between hip guides)	220mm-330mm (8.6"-12.9")	280mm-390mm (11"-15.3")
Footplate Angle	Plantarflexion 10° Dorsiflexion 10°	Plantarflexion 10° Dorsiflexion 10°
Max. Knee Contracture Accommodated	25°	25°
Tray height range (from top of hip pad)	0mm-150mm (0"-5.9")	0mm-150mm (0"-5.9")
Tray angle adjustment	45° - Prone 45° - Supine	45° - Prone 45° - Supine