



# Leckey MyWay

<b>Name</b>	
<b>Date:</b>	
<b>Email</b>	
<b>Contact Number</b>	
<b>Invoice Address</b>	
<b>Delivery Address</b>	

Quote       Order

Purchase Order	
----------------	--



	Item Code	Description	Quantity	Price
<b>Frame</b>				
<input type="checkbox"/>	173-1600	MyWay S1 MyWay Frame		XXX
<input type="checkbox"/>	173-2600	MyWay S2 MyWay Frame		XXX
<input type="checkbox"/>	173-3600	MyWay S3 MyWay Frame		XXX
<b>Harness</b>				
<input type="checkbox"/>	173-4100-04	MyWay S1 Purple MyWay Harness		XXX
<input type="checkbox"/>	173-4200-07	MyWay S2 Green MyWay Harness		XXX
<input type="checkbox"/>	173-4300-02	MyWay S3 Blue MyWay Harness		XXX
<input type="checkbox"/>	173-4400-09	MyWay S4 Red MyWay Harness		XXX
<input type="checkbox"/>	173-4500-06	MyWay S5 Orange MyWay Harness		XXX
<input type="checkbox"/>	173-4600-08	MyWay S6 Grey MyWay Harness		XXX
<input type="checkbox"/>	173-1800	MyWay Hoist Straps (one size)		XXX
<b>User Handles</b>				
<input type="checkbox"/>	173-1667	Size 1 User Handles (Pair)		XXX
<input type="checkbox"/>	173-2667	Size 2 User Handles (Pair)		XXX
<input type="checkbox"/>	173-3667	Size 3 User Handles (Pair)		XXX
<b>Adult Handles</b>				
<input type="checkbox"/>	173-1691	MyWay S1 Adult Handles		XXX
<input type="checkbox"/>	173-2691	MyWay S2 Adult Handles		XXX
<input type="checkbox"/>	173-3691	MyWay S3 Adult Handles		XXX
<input type="checkbox"/>	173-1693	MWAY Adult Connection Belt (one size) (Adult Handles required)		XXX
<b>Ankle Guides</b>				
<input type="checkbox"/>	173-1700	MyWay S1 Ankle Guides (L&R)		XXX
<input type="checkbox"/>	173-2700	MyWay S2 Ankle Guides (L&R)		XXX
<input type="checkbox"/>	173-3700	MyWay S3 Ankle Guides (L&R)		XXX
<b>Saddle</b>				
<input type="checkbox"/>	173-1671	MyWay S1 Saddle (recommended for MyWay S1 & S2 frames)		XXX
<input type="checkbox"/>	173-2671	MyWay S2 Saddle (recommended for MyWay S2 & S3 frames)		XXX
<b>MyWay Pedal (Size 2 and 3 MyWay only)</b>				
<input type="checkbox"/>	173-1750	Smart Pedal Walker with Smart Tech (one size) requires sandals		XXX
<input type="checkbox"/>	152-1600	Sandals Size 1		XXX
<input type="checkbox"/>	152-2600	Sandals Size 2		XXX
<input type="checkbox"/>	152-3600	Sandals Size 3		XXX
<input type="checkbox"/>	173-2305	Small Ankle Wraps (use with sandals)		XXX
<input type="checkbox"/>	173-3305	Large Ankle Wraps (use with sandals)		XXX

Frame Size	1	2	3
Age [approx.]	1-5	4-11	9-16
User Shoulder Height (to floor) [mm]	635-920	870-1190	1070-1490
User Saddle Height (inside leg) [mm]	220-440	400-610	550-810
Max. User Weight [kg]	25	50	80
Frame Weight [kg]	12	6.5	9.5
Upper Support Unit Weight [kg]	Size 1 frame and upper support unit do not separate	10.5	11
Frame Width	570	670	720
Frame Length	720	850	1100

Harness Size	Age [approx.]	Waist at Navel Level [cm]	Max. User Weight [kg]	Colour
1	1-4	46-56	50	Purple
2	2-6	50-60	50	Green
3	5-9	55-65	50	Blue
4	8-12	60-70	80	Red
5	12-15	70-80	80	Orange
6	Up to 16	80-90	80	Grey