



Leckey Squiggles TT

Name	
Date:	
Email	
Contact Number	
Invoice Address	
Delivery Address	

Quote Order

Purchase Order	
----------------	--



	Item Code	Description	Quantity	Price
Frame				
<input type="checkbox"/>	184-764	Squiggle + Pivot Chassis (gas operated - can be adjusted with child in frame)		XXX
<input type="checkbox"/>	184-600	Squiggle + Stander Support Shell (central support includes hip & chest support with sternum)		XXX
Cover Pack				
<input type="checkbox"/>	184-1800-09VIN	Squiggle + Stander Cover Pack (Size 1 Knees)- Red Vinyl Cover set includes 2 x support harnesses		XXX
<input type="checkbox"/>	184-1800-09	Squiggle + Stander Cover Pack (Size 1 Knees) - Red Cover set includes 2 x support harnesses		XXX
<input type="checkbox"/>	184-1800-04	Squiggle + Stander Cover Pack (Size 1 Knees)- Purple Cover set includes 2 x support harnesses		XXX
<input type="checkbox"/>	184-1800-02	Squiggle + Stander Cover Pack (Size 1 Knees)- Blue Cover set includes 2 x support harnesses		XXX
<input type="checkbox"/>	184-2800-09VIN	Squiggle + Stander Cover Pack (Size 2 Knees) - Red Vinyl Cover set includes 2 x support harnesses		XXX
<input type="checkbox"/>	184-2800-09	Squiggle +Stander Cover Pack (Size 2 Knees)- Red Cover set includes 2 x support harnesses		XXX
<input type="checkbox"/>	184-2800-04	Squiggle + Stander Cover Pack (Size 2 Knees)- Purple Cover set includes 2 x support harnesses		XXX
<input type="checkbox"/>	184-2800-02	Squiggle + Stander Cover Pack (Size 2 Knees)- Blue Cover set includes 2 x support harnesses		XXX
Tray				
<input type="checkbox"/>	118-619	Squiggles Stander Activity Tray		XXX
<input type="checkbox"/>	188-1854-05	Squiggles Black Padded Tray Cover		XXX
Foot Supports				
<input type="checkbox"/>	152-1600	Sandals Size 1		XXX
<input type="checkbox"/>	AH01	Ankle Huggers		XXX
<input type="checkbox"/>	173-2305	Small Ankle Wraps (use with sandals)		XXX
<input type="checkbox"/>	120-816	1Inch Footplate Raiser for Single Footplate		XXX
Head Supports				
<input type="checkbox"/>	137-699-05	Contoured headrest Assembly and Cushion		XXX
<input type="checkbox"/>	137-699-05VIN	Contoured headrest Assembly and Cushion - Vinyl		XXX
<input type="checkbox"/>	117-699-05	Contoured Occipital Headrest		XXX
<input type="checkbox"/>	117-699-05VIN	Contoured Occipital Headrest - Vinyl		XXX
<input type="checkbox"/>	117-851-05	Flat Headrest Assembly and Cushion		XXX
<input type="checkbox"/>	117-851-05VIN	Flat Headrest Assembly and Cushion - Vinyl		XXX
<input type="checkbox"/>	120-893-05	Flat Headrest Laterals with covers		XXX
<input type="checkbox"/>	120-893-05VIN	Flat Headrest Laterals with covers Vinyl		XXX
Targetted Training Conversion Kit				
<input type="checkbox"/>	188-680-02	Squiggle + to Squiggle TT conversion kit blue		XXX
<input type="checkbox"/>	188-680-04	Squiggle + to Squiggle TT conversion kit purple		XXX
<input type="checkbox"/>	188-680-09	Squiggle + to Squiggle TT conversion kit red		XXX
<input type="checkbox"/>	188-680-09VIN	Squiggle + to Squiggle TT conversion kit red vinyl		XXX

Age (approx)	1–5
Max User Weight	22kg/48.4lbs
User Height	Min 730mm / 28.5 inches Max 1110mm / 43.5 inches
Chest Support	Angle adjustment +/- 30° Depth adjustment 50mm / 2 inches
Chest Support Height (footplate to armpit)	Min 460mm / 18.1 inches Max 820mm / 32 inches
Chest Width (distance between laterals)	Min 160mm / 6.3 inches Max 230mm / 9 inches
Hip Support (height from footplate to mid buttock)	Min 300mm / 11.8 inches Max 530mm / 20.8 inches
Distance between midline of knees	Min 140mm / 5.5 inches Max 215mm / 8.4 inches
Footplate Angle	Plantarflexion 10° / Dorsiflexion 10°