



# Prime Engineering Kidwalk

<b>Name</b>	
<b>Date:</b>	
<b>Email</b>	
<b>Contact Number</b>	
<b>Invoice Address</b>	
<b>Delivery Address</b>	

Quote  Order

Purchase Order		
----------------	--	--



	Item Code	Description	Quantity	Price
<b>Frame</b>				
<input type="checkbox"/>	KW01	Kidwalk Size 1 (includes pelvic support swivel seat prone angle adjustment 20" quick release		XXX
<input type="checkbox"/>	KW02	Kidwalk Size 2 (includes pelvic support swivel seat prone angle adjustment 20" quick release		XXX
<b>Headrest</b>				
<input type="checkbox"/>	KW09	Kidwalk Headpad		XXX
<input type="checkbox"/>	KW10	Kidwalk Med. Contour Headrest		XXX
<b>Chest Support</b>				
<input type="checkbox"/>	KW12	Kidwalk Lateral Chest System for size 1		XXX
<input type="checkbox"/>	KW13	Kidwalk Lateral Chest System for size 2		XXX
<input type="checkbox"/>	KW20	Kidwalk Upper Body Wrap		XXX
<b>Ankle Guide</b>				
<input type="checkbox"/>	KW52	Kidwalk ankle guides Size 1 (inc attachment)		XXX
<input type="checkbox"/>	KW53	Kidwalk ankle guides Size 2 (inc attachment)		XXX
<b>Accessories</b>				
<input type="checkbox"/>	KW85	Kidwalk Directional Locks		XXX
<input type="checkbox"/>	KW03	Kidwalk Anti Rollback Device (pair)		XXX
<b>Replacement Seat</b>				
<input type="checkbox"/>	KW21	Kidwalk Wide Seat		XXX
<input type="checkbox"/>	KW07	Kidwalk 6" Extended Seat Size 1		XXX
<input type="checkbox"/>	KW23	Kidwalk 6" Extended Seat Size 2		XXX

	Size 1	Size 2
Age	1-7	5-12
Max. User Weight	29kg [65lbs]	40kg [90lbs]
Seat Height [inseam]	279mm-558mm [10.9-21.9]	381mm-711mm [15-27.9]
Hip Width	140-235mm [5.5-9.2"]	180-290mm [7-11.4"]
Chest Width	165-255mm [6.4-10"]	220-370mm [8.6-14.5"]
Thoracic Height	245-410mm [9.6-16.1"]	300-560mm [11.8-22"]