



# Leckey Everyday Activity Seat

Name	
Date:	
Email	
Contact Number	
Invoice Address	
Delivery Address	

Quote     Order

Purchase Order	
----------------	--



	Item Code	Description	Quantity	Price
<b>Seat (Black Metal)</b>				
<input type="checkbox"/>	151-6102	Size 1 seat base, backrest, and hi-low chassis with gas spring and 100mm castors (Black)		XXX
<input type="checkbox"/>	151-6204	Size 2 seat base, backrest and high-low chassis powered 100mm castors -(Black)		XXX
<input type="checkbox"/>	151-2788	Size 2 seat base, backrest and hydraulic hi low chassis with gas spring and 100mm castors (Black)		XXX
<input type="checkbox"/>	151-6202	Size 2 seat base, backrest and hi-low chassis with gas spring and 100mm castors (Black)		XXX
<input type="checkbox"/>	151-7204	Size 2 Extended seat base, backrest and high-low chassis powered 100mm castors -Black		XXX
<input type="checkbox"/>	151-7202	Size 2 Extended seat base, backrest and hi-low chassis with gas spring and 100mm castors (Black)		XXX
<input type="checkbox"/>	151-6340	Size 3 seat base, backrest and a Hi Low Powered Chassis (includes pushbar) - Black		XXX
<input type="checkbox"/>	151-6302	Size 3 seat base, backrest and a hydraulic chassis with 125mm castors -Black		XXX
<b>Seat (White Metal)</b>				
<input type="checkbox"/>	151-6103	Size 1 seat base, backrest, and hi-low chassis with gas spring and 100mm castors (White)		XXX
<input type="checkbox"/>	151-6203	Size 2 seat base, backrest and hi-low chassis with gas spring and 100mm castors (White)		XXX
<input type="checkbox"/>	151-7203	Size 2 Extended seat base, backrest and hi-low chassis with gas spring and 100mm castors (White)		XXX
<input type="checkbox"/>	151-7205	Size 2 Extended seat base, backrest and high-low chassis powered 100mm castors -White		XXX
<input type="checkbox"/>	151-6303	Size 3 seat base, backrest and a hydraulic chassis with 125mm castors -White		XXX
<b>Covers</b>				
<input type="checkbox"/>	151-6104-07VIN	Size 1 cushions and covers - black vinyl		XXX
<input type="checkbox"/>	151-6104-05	Size 1 cushions and covers - ebony grey		XXX
<input type="checkbox"/>	151-6104-03VIN	Size 1 cushions and covers - blue vinyl		XXX
<input type="checkbox"/>	151-6104-03	Size 1 cushions and covers - blue		XXX
<input type="checkbox"/>	151-6104-01	Size 1 cushions and covers - green		XXX
<input type="checkbox"/>	151-6206-07VIN	Size 2 cushions and covers - Black vinyl		XXX
<input type="checkbox"/>	151-6206-05	Size 2 cushions and covers - Ebony Grey		XXX
<input type="checkbox"/>	151-6206-03VIN	Size 2 cushions and covers - Blue vinyl		XXX
<input type="checkbox"/>	151-6206-03	Size 2 cushions and covers - Blue		XXX
<input type="checkbox"/>	151-6206-01	Size 2 cushions and covers - Green		XXX
<input type="checkbox"/>	151-6240-07VIN	Size 2 Extended cushions and covers - Black Vinyl (with sacral foam spacer pads)		XXX
<input type="checkbox"/>	151-6240-05	Size 2 Extended cushions and covers - Ebony Grey base (with sacral foam spacer pads)		XXX
<input type="checkbox"/>	151-6240-03VIN	Size 2 Extended cushions and covers - Blue base (with sacral foam spacer pads)		XXX
<input type="checkbox"/>	151-6240-03	Size 2 Extended cushions and covers - Blue base (with sacral foam spacer pads)		XXX
<input type="checkbox"/>	151-6240-01	Size 2 Extended cushions and covers - Green base (with sacral foam spacer pads)		XXX
<input type="checkbox"/>	151-6304-07VIN	Size 3 cushions and covers - Black vinyl		XXX
<input type="checkbox"/>	151-6304-05	Size 3 cushions and covers - Ebony Grey		XXX
<input type="checkbox"/>	151-6304-03VIN	Size 3 cushions and covers - Blue vinyl		XXX
<input type="checkbox"/>	151-6304-03	Size 3 cushions and covers - Blue		XXX
<input type="checkbox"/>	151-6304-01	Size 3 cushions and covers - Green		XXX
<b>Pommel (Black Metal)</b>				
<input type="checkbox"/>	151-6120-07VIN	Wedge shape pommel S1 - Black Vinyl		XXX
<input type="checkbox"/>	151-6120-05	Wedge shape pommel S1 - Ebony Grey		XXX
<input type="checkbox"/>	151-6118-07VIN	Size 1 Disc Pommel - Black Vinyl		XXX
<input type="checkbox"/>	151-6118-05	Size 1 Disc Pommel - Ebony Grey		XXX
<input type="checkbox"/>	151-6229-07VIN	Wedge shape pommel S2 - Black Vinyl		XXX
<input type="checkbox"/>	151-6229-05	Wedge shape pommel S2 - Ebony Grey		XXX
<input type="checkbox"/>	151-6227-07VIN	Size 2 Disc Pommel - Black Vinyl		XXX
<input type="checkbox"/>	151-6227-05	Size 2 Disc Pommel - Ebony Grey		XXX
<input type="checkbox"/>	151-6334-07VIN	Wedge shape pommel S3 - Black Vinyl		XXX
<input type="checkbox"/>	151-6334-05	Wedge shape pommel S3 - Ebony Grey		XXX
<input type="checkbox"/>	151-6322-07VIN	Size 3 Disc Pommel - Black Vinyl		XXX
<input type="checkbox"/>	151-6322-05	Size 3 Disc Pommel - Ebony Grey		XXX
<b>Pommel (White Metal)</b>				

<input type="checkbox"/>	151-6119-05	Size 1 Disc Pommel - Ebony Grey (White)		XXX
<input type="checkbox"/>	151-6228-05	Size 2 Disc Pommel - Ebony Grey (White)		XXX
<input type="checkbox"/>	151-6323-05	Size 3 Disc Pommel - Ebony Grey (White)		XXX
<input type="checkbox"/>	151-6121-05	Wedge shape pommel S1 - Ebony Grey (White)		XXX
<input type="checkbox"/>	151-6230-05	Wedge shape pommel S2 - Ebony Grey (White)		XXX
<input type="checkbox"/>	151-6335-05	Wedge shape pommel S3 - Ebony Grey (White)		XXX
<b>Armrest (Black Metal)</b>				
<input type="checkbox"/>	151-6117	Size 1 Armrest (Pair)		XXX
<input type="checkbox"/>	151-6226	Size 2 Armrest (Pair)		XXX
<input type="checkbox"/>	151-3747	Size 3 armrests pair		XXX
<b>Armrest (White Metal)</b>				
<input type="checkbox"/>	151-6124	Size 1 Armrests (White)		XXX
<input type="checkbox"/>	151-6234	Size 2 Armrest (Pair) White		XXX
<input type="checkbox"/>	151-6328	Size 3 Armrest (Pair) White		XXX
<b>Footplate</b>				
<input type="checkbox"/>	151-1745	Size 1 One piece footplate (Black)		XXX
<input type="checkbox"/>	151-2745	Size 2 One piece footplate (Black)		XXX
<input type="checkbox"/>	151-3631	Size 3 One Piece Padded Footplate (Black)		XXX
<input type="checkbox"/>	151-3745	Size 3 One piece footplate (Black)		XXX
<input type="checkbox"/>	151-6123	S1 One piece footplate (White)		XXX
<input type="checkbox"/>	151-6223	S2 One piece footplate (White)		XXX
<input type="checkbox"/>	151-6327	S3 One piece footplate (White)		XXX
<input type="checkbox"/>	151-3783	Lower Leg Strap		XXX
<b>Tray</b>				
<input type="checkbox"/>	151-1644	Size 1 Tray Plastic		XXX
<input type="checkbox"/>	151-2644	Size 2 Tray Plastic		XXX
<input type="checkbox"/>	151-3767	Size 3 Tray assembly black (additional 20mm seat width)		XXX
<input type="checkbox"/>	151-6122	S1 Tray Wooden flipaway		XXX
<input type="checkbox"/>	151-6231	S2 Tray Wooden flipaway		XXX
<input type="checkbox"/>	151-6326	S3 Tray Wooden Flipaway		XXX
<b>Laterals (Black Metal)</b>				
<input type="checkbox"/>	151-6108-07VIN	Size 1 flip away PU chest laterals pair - black vinyl		XXX
<input type="checkbox"/>	151-6108-05	Size 1 flip away PU chest laterals pair - ebony grey		XXX
<input type="checkbox"/>	151-6106-07VIN	Size 1 Rigid PU chest laterals pair - black vinyl		XXX
<input type="checkbox"/>	151-6106-05	Size 1 Rigid PU chest laterals pair - Ebony Grey		XXX
<input type="checkbox"/>	151-1757	Size 1 complex flipaway laterals pair		XXX
<input type="checkbox"/>	151-6212-07VIN	S2 Flip away PU chest laterals pair - black vinyl		XXX
<input type="checkbox"/>	151-6212-05	S2 Flip away PU chest laterals pair - ebony grey		XXX
<input type="checkbox"/>	151-6210-07VIN	S2 Rigid PU chest laterals pair - black vinyl		XXX
<input type="checkbox"/>	151-6210-05	S2 Rigid PU chest laterals pair - ebony Grey		XXX
<input type="checkbox"/>	151-2757	Size 2 complex flipaway laterals pair		XXX
<input type="checkbox"/>	151-6308-07VIN	S3 Flip away PU chest laterals pair - Black Vinyl		XXX
<input type="checkbox"/>	151-6308-05	S3 Flip away PU chest laterals pair - ebony grey		XXX
<input type="checkbox"/>	151-6306-07VIN	S3 Rigid PU chest laterals pair - black vinyl		XXX
<input type="checkbox"/>	151-6306-05	S3 Rigid PU chest laterals pair - ebony grey		XXX
<input type="checkbox"/>	151-3757	Size 3 complex flipaway laterals		XXX
<b>Laterals (White Metal)</b>				
<input type="checkbox"/>	151-6107-05	Size 1 Rigid PU chest laterals pair - ebony grey (White Metal)		XXX
<input type="checkbox"/>	151-6211-05	S2 Rigid PU chest laterals pair - ebony grey (White Metal)		XXX
<input type="checkbox"/>	151-6307-05	S3 Rigid PU chest laterals pair - ebony grey (White Metal)		XXX
<input type="checkbox"/>	151-6109-05	Size 1 flip away PU chest laterals pair - ebony grey (White Metal)		XXX
<input type="checkbox"/>	151-6213-05	S2 Flip away PU chest laterals pair - ebony grey (White Metal)		XXX
<input type="checkbox"/>	151-6309-05	S3 Flip away PU chest laterals pair - ebony grey (White Metal)		XXX
<b>Everyday Seat Headrest</b>				
<input type="checkbox"/>	151-6209-05VIN	Size 1 / 2 Flat headrest laterals - Ebony Grey Vinyl (only available with 01VIN, 07VIN & 04VIN c		XXX
<input type="checkbox"/>	151-6209-05	S 1 / 2 Flat headrest laterals - Ebony Grey (only available with 01,03 & 04 cushion packs)		XXX
<input type="checkbox"/>	151-6208-05VIN	Size 1 & 2 Flat headrest - Ebony Grey Vinyl (only available with 01VIN, 07VIN & 04VIN cushion		XXX
<input type="checkbox"/>	151-6208-05	Size 1 & 2 Flat headrest - Ebony Grey (only available with 01,03 & 04 cushion packs)		XXX
<input type="checkbox"/>	151-6105-05VIN	Size 1 Contoured Headrest and Cover - Grey Vinyl (only available with 01, 07 & 04 cushion pac		XXX
<input type="checkbox"/>	151-6105-05	Size 1 Contoured Headrest and Cover - Ebony Grey (only available with 01, 03 & 04 cushion pa		XXX
<input type="checkbox"/>	151-6207-07VIN	Size 2 Contoured Headrest and Cover - Black Vinyl (only available with -03VIN cushion pack)		XXX
<input type="checkbox"/>	151-6207-05	Size 2 Contoured Headrest and Cover - Ebony Grey (only available with -05 cushion pack)		XXX
<input type="checkbox"/>	151-3751-05	Size 3 Contoured Headrest - (black neoprene)		XXX
<b>Other Headrests (Requires Bracket)</b>				
<input type="checkbox"/>	151-1756	Universal headrest attachment bracket		XXX
<input type="checkbox"/>	137-699-05	Contoured headrest Assembly and Cushion		XXX
<input type="checkbox"/>	137-699-05VIN	Contoured headrest Assembly and Cushion - Vinyl		XXX
<input type="checkbox"/>	117-699-05	Contoured Occipital Headrest Foam & Cover		XXX
<input type="checkbox"/>	117-699-05VIN	Contoured Occipital Headrest Foam & Cover - Vinyl		XXX

<input type="checkbox"/>	117-851-05	Flat headrest Assembly and Cushion		XXX
<input type="checkbox"/>	117-851-05VIN	Flat headrest Assembly and Cushion - Vinyl		XXX
<input type="checkbox"/>	120-893-05	Flat Headrest Laterals with covers		XXX
<input type="checkbox"/>	120-893-05VIN	Flat Headrest Laterals with covers Vinyl		XXX
<b>Side Pads</b>				
<input type="checkbox"/>	151-6115-07VIN	Size 1 side pads pair - Black Vinyl		XXX
<input type="checkbox"/>	151-6115-05	Size 1 side pads pair - Ebony Grey		XXX
<input type="checkbox"/>	151-6222-07VIN	Size 2 side pads pair - Black Vinyl		XXX
<input type="checkbox"/>	151-6222-05	Size 2 side pads pair - Ebony Grey		XXX
<input type="checkbox"/>	151-6317-07VIN	Size 3 side pads pair - Black Vinyl		XXX
<input type="checkbox"/>	151-6317-05	Size 3 side pads pair - Ebony Grey		XXX
<input type="checkbox"/>	151-6116-05	Size 1 side pads pair - Ebony Grey (White Metal)		XXX
<input type="checkbox"/>	151-6223-05	Size 2 side pads pair - Ebony Grey (White Metal)		XXX
<input type="checkbox"/>	151-6318-05	Size 3 side pads pair - Ebony Grey (White Metal)		XXX
<b>Protraction Pad</b>				
<input type="checkbox"/>	151-6218-05VIN	Size 1/2 Protraction pad pair Ebony Grey Vinyl		XXX
<input type="checkbox"/>	151-6218-05	Size 1/2 Protraction pad pair - Ebony Grey		XXX
<input type="checkbox"/>	151-6314-05VIN	Size 3 Protraction Pad pair - Ebony Grey Vinyl		XXX
<input type="checkbox"/>	151-6314-05	Size 3 Protraction Pad pair - Ebony Grey		XXX
<input type="checkbox"/>	151-6219-05	Size 1/2 Protraction pad pair Ebony Grey Vinyl (White Metal)		XXX
<input type="checkbox"/>	151-6315-05	Size 3 Protraction Pad pair - Ebony Grey (White Metal)		XXX
<b>Abduction Pad</b>				
<input type="checkbox"/>	151-6224-07VIN	Size 2 Abduction Side Pad - Black Vinyl		XXX
<input type="checkbox"/>	151-6224-05	Size 2 Abduction Side Pad - Ebony Grey		XXX
<input type="checkbox"/>	151-6320-07VIN	Size 3 Abduction Side Pad - Black Vinyl		XXX
<input type="checkbox"/>	151-6320-05	Size 3 Abduction Side Pad - Ebony Grey		XXX
<input type="checkbox"/>	151-6225-05	Size 2 Abduction Side Pad - Ebony Grey (White Metal)		XXX
<input type="checkbox"/>	151-6321-05	Size 3 Abduction Side Pad - Ebony Grey (White Metal)		XXX
<b>Chest Support</b>				
<input type="checkbox"/>	137-808-05	Oval Chest Harness - Black		XXX
<input type="checkbox"/>	137-808-05VIN	Oval Chest Harness - Black Vinyl		XXX
<input type="checkbox"/>	155-N001-L600	Small Butterfly Harness & attachment straps		XXX
<input type="checkbox"/>	155-N001-L601	Medium Butterfly Harness & attachment straps		XXX
<input type="checkbox"/>	155-N001-L602	Large Butterfly Harness & attachment straps		XXX
<input type="checkbox"/>	155-N001-L603	X-Large Butterfly Harness & attachment straps		XXX
<input type="checkbox"/>	133-625	EDS size 3 Chest Pad		XXX
<b>Foot Support</b>				
<input type="checkbox"/>	152-1600	Sandals Size 1		XXX
<input type="checkbox"/>	152-2600	Sandals Size 2		XXX
<input type="checkbox"/>	152-3600	Sandals Size 3		XXX
<input type="checkbox"/>	AH01	Ankle Huggers		XXX
<input type="checkbox"/>	173-2305	Small Ankle Wraps (use with sandals)		XXX
<input type="checkbox"/>	173-3305	Large Ankle Wraps (use with sandals)		XXX
<input type="checkbox"/>	117-787	1 Inch Sandal Raisers (pair) for split footplate		XXX
<b>Pelvic Support</b>				
<input type="checkbox"/>	151-1660	4 point pelvic belt - Size 1		XXX
<input type="checkbox"/>	151-2660	4 point pelvic belt - Size 2		XXX
<input type="checkbox"/>	151-3660	4 point pelvic belt - Size 3		XXX
<input type="checkbox"/>	117-1860	2-point belt - Size 1		XXX
<input type="checkbox"/>	117-2860	2-point belt - Size 2		XXX
<input type="checkbox"/>	117-3860	2-point belt - Size 3		XXX

Technical Details			
	Size 1	Size 2	Size 3
<b>Maximum load</b>	30kgs 66lbs	50kg 110lbs	100kgs 220lbs
<b>Weight - seat and hi-low chassis</b>	29.5kg 65lbs	34kg 75lbs	49kg 108lbs
<b>Seat height – hi-low chassis</b>	250-600 9.8-23.6	320-600 12.6-23.6	450-600 17.7-23.6
<b>Seat width</b>	160-280 6.3-11	260-380 10.2-15	340-500 13.4-19.7
<b>Seat depth (standard)</b>	150-280 5.9-11	230-380 9-15	330-480 13-1109
<b>Seat depth (extended)</b>	150-280 5.9-11	300-450 11.8-17.7	330-480 13-1109
<b>Lateral chest distance</b>	150-220 5.9-10.7	170-240 6.7-9.4	240-370 9.4-14.6
<b>Back height</b>	300-400 11.8-15.7	400-520 15.7-20.5	450-630 17.7-24.8
<b>Back angle</b>	30° supine 10° prone	30° supine 10° prone	30° supine 10° prone
<b>Armrest height</b>	150-250 5.9-9.8	200-300 7.9-11.8	180-325 7-12.7
<b>Lower leg length</b>	150-250 5.9-9.8	250-450 9.8-17.7	250-450 9.8-17.7

<b>Lower leg angle</b>	-10° to +50°	-10° to +50°	-10° to +50°
<b>Footplate angle</b>	-15° to +10°	-15° to +10°	-15° to +10°
<b>Tilt-in-space</b>	30° supine 10° prone	30° supine 10° prone	30° supine 10° prone