



Firefly WeGo CE

Name	
Date:	
Email	
Contact Number	
Invoice Address	
Delivery Address	

Quote Order

Purchase Order	
----------------	--



	Item Code	Description	Quantity	Price
<input type="checkbox"/>	182-500	WEGO WeGo Push Chair (includes all accessories)		XXX
<input type="checkbox"/>	182-501	WEGO WeGo Accessories		XXX
Firefly GoTo Seat - Individual				
<input type="checkbox"/>	167-501-09-S	GoTo Size 1 - Standard Headrest - Red Vinyl		XXX
<input type="checkbox"/>	167-501-09-C	GoTo Size 1 - Advanced Headrest - Red Vinyl		XXX
<input type="checkbox"/>	167-501-04-C	GoTo Size 1 - Advanced Headrest - Pink / Purple		XXX
<input type="checkbox"/>	167-501-04-S	GoTo Size 1 - Standard Headrest - Pink / Purple		XXX
<input type="checkbox"/>	167-501-05-C	GoTo Size 1 - Advanced Headrest - Light Blue / Dark Blue		XXX
<input type="checkbox"/>	167-501-05-S	GoTo Size 1 - Standard Headrest - Light Blue / Dark Blue		XXX
<input type="checkbox"/>	167-502-09-S	GoTo Size 2 - Standard Headrest - Red Vinyl		XXX
<input type="checkbox"/>	167-502-09-C	GoTo Size 2 - Advanced Headrest - Red Vinyl		XXX
<input type="checkbox"/>	167-502-05-S	GoTo Size 2 - Standard Headrest - Light Blue / Dark Blue		XXX
<input type="checkbox"/>	167-502-05-C	GoTo Size 2 - Advanced Headrest - Light Blue / Dark Blue		XXX
<input type="checkbox"/>	167-502-04-S	GoTo Size 2 - Standard Headrest - Pink / Purple		XXX
<input type="checkbox"/>	167-502-04-C	GoTo Size 2 - Advanced Headrest - Pink / Purple		XXX
<input type="checkbox"/>	167-101	Size 1 GoTo Cushion & Roll		XXX
<input type="checkbox"/>	167-201	Size 2 GoTo Cushion & Roll		XXX
<input type="checkbox"/>	167-101-09	Size 1 GoTo Cushion & Roll - Red Vinyl		XXX
<input type="checkbox"/>	167-201-09	Size 1 GoTo Cushion & Roll - Red Vinyl		XXX
<input type="checkbox"/>	176-500	Floorsitter		XXX
Firefly GoTo Seat - Package				
<input type="checkbox"/>	167-501-09-SF-CU	GoTo Size 1 - Std Headrest - Red Vinyl with Floorsitter with Pressure Cushion		XXX
<input type="checkbox"/>	167-501-09-CF-CU	GoTo Size 1 - Advanced Headrest - Red Vinyl with Floorsitter with Pressure Cushion		XXX
<input type="checkbox"/>	167-501-05-SF-CU	GoTo Size 1 - Std Headrest - Light Blue / Dark Blue with Floorsitter with Pressure Cushion		XXX
<input type="checkbox"/>	167-501-05-CF-CU	GoTo Size 1 - Advanced Headrest - Light Blue / Dark Blue with Floorsitter with Pressure Cushion		XXX
<input type="checkbox"/>	167-501-04-SF-CU	GoTo Size 1 - Std Headrest - Pink / Purple with Floorsitter with Pressure Cushion		XXX
<input type="checkbox"/>	167-501-04-CF-CU	GoTo Size 1 - Advanced Headrest - Pink / Purple with Floorsitter with Pressure Cushion		XXX
<input type="checkbox"/>	167-502-09-SF-CU	GoTo Size 2 - Std Headrest - Red Vinyl with Floorsitter with Pressure Cushion		XXX
<input type="checkbox"/>	167-502-09-CF-CU	GoTo Size 2 - Advanced Headrest - Red Vinyl with Floorsitter with Pressure Cushion		XXX
<input type="checkbox"/>	167-502-05-SF-CU	GoTo Size 2 - Std Headrest - Light Blue / Dark Blue with Floorsitter with Pressure Cushion		XXX
<input type="checkbox"/>	167-502-05-CF-CU	GoTo Size 2 - Advanced Headrest - Light Blue / Dark Blue with Floorsitter with Pressure Cushion		XXX
<input type="checkbox"/>	167-502-04-SF-CU	GoTo Size 2 - Std Headrest - Pink / Purple with Floorsitter with Pressure Cushion		XXX
<input type="checkbox"/>	167-502-04-CF-CU	GoTo Size 2 - Advanced Headrest - Pink / Purple with Floorsitter with Pressure Cushion		XXX

GoTo Seat	Size 1	Size 2
Base to Top of Headrest	380mm-540mm [14"-21"]	510mm-680mm [20"-27"]
Lateral Support Width	190mm-260mm [7.5"-10"]	210mm-300mm [8"-12"]
Seat Depth	175mm [7"]	244mm [9.5"]
Max. User Weight	15kg [33lbs]	30kg [66lbs]
Age Range	1-3 years	3-8 years

WeGo	
Footprint Length (Collapsed)	116cm
Footprint Height (Collapsed)	38cm
WeGo Height	112cm
Seat Depth	26cm
Seat Height	50cm
Seat Width	35cm
Footrest Shortest setting	25cm
Footrest Longest setting	33cm
WeGo Weight	11kg
WeGo Weight with GoTo	13kg