



Leckey Everyday Activity Seat



CE

Name	
Date:	
Email	
Contact Number	
Invoice Address	
Delivery Address	


Quote Order

Purchase Order	
----------------	--



	Item Code	Description	Quantity	Price
Seat (Black Metal)				
				
<input type="checkbox"/>	151-6102	Size 1 seat base, backrest, and hi-low chassis with gas spring and 100mm castors (Black)		XXX
<input type="checkbox"/>	151-6204	Size 2 seat base, backrest and high-low chassis powered 100mm castors -(Black)		XXX
<input type="checkbox"/>	151-2788	Size 2 seat base, backrest and hydraulic hi low chassis with gas spring and 100mm castors (Black)		XXX
<input type="checkbox"/>	151-6202	Size 2 seat base, backrest and hi-low chassis with gas spring and 100mm castors (Black)		XXX
<input type="checkbox"/>	151-7204	Size 2 Extended seat base, backrest and high-low chassis powered 100mm castors -Black		XXX
<input type="checkbox"/>	151-7202	Size 2 Extended seat base, backrest and hi-low chassis with gas spring and 100mm castors (Black)		XXX
<input type="checkbox"/>	151-6340	Size 3 seat base, backrest and a Hi Low Powered Chassis (includes pushbar) - Black		XXX
<input type="checkbox"/>	151-6302	Size 3 seat base, backrest and a hydraulic chassis with 125mm castors -Black		XXX
Seat (White Metal)				
				
<input type="checkbox"/>	151-6103	Size 1 seat base, backrest, and hi-low chassis with gas spring and 100mm castors (White)		XXX
<input type="checkbox"/>	151-6203	Size 2 seat base, backrest and hi-low chassis with gas spring and 100mm castors (White)		XXX
<input type="checkbox"/>	151-7203	Size 2 Extended seat base, backrest and hi-low chassis with gas spring and 100mm castors (White)		XXX
<input type="checkbox"/>	151-7205	Size 2 Extended seat base, backrest and high-low chassis powered 100mm castors -White		XXX
<input type="checkbox"/>	151-6303	Size 3 seat base, backrest and a hydraulic chassis with 125mm castors -White		XXX

Covers




<input type="checkbox"/>	151-6104-07VIN	Size 1 cushions and covers - black vinyl	XXX
<input type="checkbox"/>	151-6104-05	Size 1 cushions and covers - ebony grey	XXX
<input type="checkbox"/>	151-6104-03VIN	Size 1 cushions and covers - blue vinyl	XXX
<input type="checkbox"/>	151-6104-03	Size 1 cushions and covers - blue	XXX
<input type="checkbox"/>	151-6104-01	Size 1 cushions and covers - green	XXX
<input type="checkbox"/>	151-6206-07VIN	Size 2 cushions and covers - Black vinyl	XXX
<input type="checkbox"/>	151-6206-05	Size 2 cushions and covers - Ebony Grey	XXX
<input type="checkbox"/>	151-6206-03VIN	Size 2 cushions and covers - Blue vinyl	XXX
<input type="checkbox"/>	151-6206-03	Size 2 cushions and covers - Blue	XXX
<input type="checkbox"/>	151-6206-01	Size 2 cushions and covers - Green	XXX
<input type="checkbox"/>	151-6240-07VIN	Size 2 Extended cushions and covers - Black Vinyl (with sacral foam spacer pads)	XXX
<input type="checkbox"/>	151-6240-05	Size 2 Extended cushions and covers - Ebony Grey base (with sacral foam spacer pads)	XXX
<input type="checkbox"/>	151-6240-03VIN	Size 2 Extended cushions and covers - Blue base (with sacral foam spacer pads)	XXX
<input type="checkbox"/>	151-6240-03	Size 2 Extended cushions and covers - Blue base (with sacral foam spacer pads)	XXX
<input type="checkbox"/>	151-6240-01	Size 2 Extended cushions and covers - Green base (with sacral foam spacer pads)	XXX
<input type="checkbox"/>	151-6304-07VIN	Size 3 cushions and covers - Black vinyl	XXX
<input type="checkbox"/>	151-6304-05	Size 3 cushions and covers - Ebony Grey	XXX
<input type="checkbox"/>	151-6304-03VIN	Size 3 cushions and covers - Blue vinyl	XXX
<input type="checkbox"/>	151-6304-03	Size 3 cushions and covers - Blue	XXX
<input type="checkbox"/>	151-6304-01	Size 3 cushions and covers - Green	XXX

Pommel (Black Metal)



Disc




Wedge

<input type="checkbox"/>	151-6120-07VIN	Wedge shape pommel S1 - Black Vinyl		XXX
<input type="checkbox"/>	151-6120-05	Wedge shape pommel S1 - Ebony Grey		XXX
<input type="checkbox"/>	151-6118-07VIN	Size 1 Disc Pommel - Black Vinyl		XXX
<input type="checkbox"/>	151-6118-05	Size 1 Disc Pommel - Ebony Grey		XXX
<input type="checkbox"/>	151-6229-07VIN	Wedge shape pommel S2 - Black Vinyl		XXX
<input type="checkbox"/>	151-6229-05	Wedge shape pommel S2 - Ebony Grey		XXX
<input type="checkbox"/>	151-6227-07VIN	Size 2 Disc Pommel - Black Vinyl		XXX
<input type="checkbox"/>	151-6227-05	Size 2 Disc Pommel - Ebony Grey		XXX
<input type="checkbox"/>	151-6334-07VIN	Wedge shape pommel S3 - Black Vinyl		XXX
<input type="checkbox"/>	151-6334-05	Wedge shape pommel S3 - Ebony Grey		XXX
<input type="checkbox"/>	151-6322-07VIN	Size 3 Disc Pommel - Black Vinyl		XXX
<input type="checkbox"/>	151-6322-05	Size 3 Disc Pommel - Ebony Grey		XXX

Pommel (White Metal)



Disc



Wedge

<input type="checkbox"/>	151-6119-05	Size 1 Disc Pommel - Ebony Grey (White)		XXX
<input type="checkbox"/>	151-6228-05	Size 2 Disc Pommel - Ebony Grey (White)		XXX
<input type="checkbox"/>	151-6323-05	Size 3 Disc Pommel - Ebony Grey (White)		XXX
<input type="checkbox"/>	151-6121-05	Wedge shape pommel S1 - Ebony Grey (White)		XXX
<input type="checkbox"/>	151-6230-05	Wedge shape pommel S2 - Ebony Grey (White)		XXX
<input type="checkbox"/>	151-6335-05	Wedge shape pommel S3 - Ebony Grey (White)		XXX

Armrest (Black Metal)



<input type="checkbox"/>	151-6117	Size 1 Armrest (Pair)		XXX
<input type="checkbox"/>	151-6226	Size 2 Armrest (Pair)		XXX
<input type="checkbox"/>	151-3747	Size 3 armrests pair		XXX

Armrest (White Metal)



<input type="checkbox"/>	151-6124	Size 1 Armrests (White)		XXX
<input type="checkbox"/>	151-6234	Size 2 Armrest (Pair) White		XXX
<input type="checkbox"/>	151-6328	Size 3 Armrest (Pair) White		XXX

Footplate



Black Size 1/2 Footplate



Black Size 3 Footplate



White Footplate

<input type="checkbox"/>	151-1745	Size 1 One piece footplate (Black)	XXX
<input type="checkbox"/>	151-2745	Size 2 One piece footplate (Black)	XXX
<input type="checkbox"/>	151-3631	Size 3 One Piece Padded Footplate (Black)	XXX
<input type="checkbox"/>	151-3745	Size 3 One piece footplate (Black)	XXX
<input type="checkbox"/>	151-6123	S1 One piece footplate (White)	XXX
<input type="checkbox"/>	151-6223	S2 One piece footplate (White)	XXX
<input type="checkbox"/>	151-6327	S3 One piece footplate (White)	XXX
<input type="checkbox"/>	151-3783	Lower Leg Strap	XXX

Tray



Grey Size 1/2 Tray



Black Size 3 Tray



Wooden Flipaway Tray

<input type="checkbox"/>	151-1644	Size 1 Tray Plastic	XXX
<input type="checkbox"/>	151-2644	Size 2 Tray Plastic	XXX
<input type="checkbox"/>	151-3767	Size 3 Tray assembly black (additional 20mm seat width)	XXX
<input type="checkbox"/>	151-6122	Size 1 Tray Wooden flipaway	XXX
<input type="checkbox"/>	151-6231	Size 2 Tray Wooden flipaway	XXX
<input type="checkbox"/>	151-6326	Size 3 Tray Wooden Flipaway	XXX

Laterals (Black Metal)



Rigid PU Chest Laterals



Flip-away PU Chest Laterals



Complex Flip-Away Chest
Laterals

<input type="checkbox"/>	151-6108-07VIN	Size 1 flip away PU chest laterals pair - black vinyl		XXX
<input type="checkbox"/>	151-6108-05	Size 1 flip away PU chest laterals pair - ebony grey		XXX
<input type="checkbox"/>	151-6106-07VIN	Size 1 Rigid PU chest laterals pair - black vinyl		XXX
<input type="checkbox"/>	151-6106-05	Size 1 Rigid PU chest laterals pair - Ebony Grey		XXX
<input type="checkbox"/>	151-1757	Size 1 complex flipaway laterals pair		XXX
<input type="checkbox"/>	151-6212-07VIN	S2 Flip away PU chest laterals pair - black vinyl		XXX
<input type="checkbox"/>	151-6212-05	S2 Flip away PU chest laterals pair - ebony grey		XXX
<input type="checkbox"/>	151-6210-07VIN	S2 Rigid PU chest laterals pair - black vinyl		XXX
<input type="checkbox"/>	151-6210-05	S2 Rigid PU chest laterals pair - ebony Grey		XXX
<input type="checkbox"/>	151-2757	Size 2 complex flipaway laterals pair		XXX
<input type="checkbox"/>	151-6308-07VIN	S3 Flip away PU chest laterals pair - Black Vinyl		XXX
<input type="checkbox"/>	151-6308-05	S3 Flip away PU chest laterals pair - ebony grey		XXX
<input type="checkbox"/>	151-6306-07VIN	S3 Rigid PU chest laterals pair - black vinyl		XXX
<input type="checkbox"/>	151-6306-05	S3 Rigid PU chest laterals pair - ebony grey		XXX
<input type="checkbox"/>	151-3757	Size 3 complex flipaway laterals		XXX

Laterals (White Metal)



Rigid PU Chest Laterals



Flip-away PU Chest Laterals

<input type="checkbox"/>	151-6107-05	Size 1 Rigid PU chest laterals pair - ebony grey (White Metal)		XXX
<input type="checkbox"/>	151-6211-05	S2 Rigid PU chest laterals pair - ebony grey (White Metal)		XXX
<input type="checkbox"/>	151-6307-05	S3 Rigid PU chest laterals pair - ebony grey (White Metal)		XXX
<input type="checkbox"/>	151-6109-05	Size 1 flip away PU chest laterals pair - ebony grey (White Metal)		XXX
<input type="checkbox"/>	151-6213-05	S2 Flip away PU chest laterals pair - ebony grey (White Metal)		XXX
<input type="checkbox"/>	151-6309-05	S3 Flip away PU chest laterals pair - ebony grey (White Metal)		XXX

Everyday Seat Headrest



Contoured Occipital



Flat



Laterals (for use with flat)



Size 3 Contoured

<input type="checkbox"/>	151-6209-05VIN	Size 1 / 2 Flat headrest laterals - Ebony Grey Vinyl (only available with 01VIN, 07VIN & 04VIN c	XXX
<input type="checkbox"/>	151-6209-05	S 1 / 2 Flat headrest laterals - Ebony Grey (only available with 01,03 & 04 cushion packs)	XXX
<input type="checkbox"/>	151-6208-05VIN	Size 1 & 2 Flat headrest - Ebony Grey Vinyl (only available with 01VIN, 07VIN & 04VIN cushion	XXX
<input type="checkbox"/>	151-6208-05	Size 1 & 2 Flat headrest - Ebony Grey (only available with 01,03 & 04 cushion packs)	XXX
<input type="checkbox"/>	151-6105-05VIN	Size 1 Contoured Headrest and Cover - Grey Vinyl (only available with 01, 07 & 04 cushion pac	XXX
<input type="checkbox"/>	151-6105-05	Size 1 Contoured Headrest and Cover - Ebony Grey (only available with 01, 03 & 04 cushion pa	XXX
<input type="checkbox"/>	151-6207-07VIN	Size 2 Contoured Headrest and Cover - Black Vinyl (only available with -03VIN cushion pack)	XXX
<input type="checkbox"/>	151-6207-05	Size 2 Contoured Headrest and Cover - Ebony Grey (only available with -05 cushion pack)	XXX
<input type="checkbox"/>	151-3751-05	Size 3 Contoured Headrest - (black neoprene)	XXX

Other Headrests (Requires Bracket)



Contoured



Contoured Occipital



Flat



Laterals (for use with flat)

<input type="checkbox"/>	151-1756	Universal headrest attachment bracket	XXX
<input type="checkbox"/>	137-699-05	Contoured headrest Assembly and Cushion	XXX
<input type="checkbox"/>	137-699-05VIN	Contoured headrest Assembly and Cushion - Vinyl	XXX
<input type="checkbox"/>	117-699-05	Contoured Occipital Headrest Foam & Cover	XXX
<input type="checkbox"/>	117-699-05VIN	Contoured Occipital Headrest Foam & Cover - Vinyl	XXX
<input type="checkbox"/>	117-851-05	Flat headrest Assembly and Cushion	XXX
<input type="checkbox"/>	117-851-05VIN	Flat headrest Assembly and Cushion - Vinyl	XXX
<input type="checkbox"/>	120-893-05	Flat Headrest Laterals with covers	XXX
<input type="checkbox"/>	120-893-05VIN	Flat Headrest Laterals with covers Vinyl	XXX

Side Pads



<input type="checkbox"/>	151-6115-07VIN	Size 1 side pads pair - Black Vinyl	XXX
<input type="checkbox"/>	151-6115-05	Size 1 side pads pair - Ebony Grey	XXX
<input type="checkbox"/>	151-6222-07VIN	Size 2 side pads pair - Black Vinyl	XXX
<input type="checkbox"/>	151-6222-05	Size 2 side pads pair - Ebony Grey	XXX
<input type="checkbox"/>	151-6317-07VIN	Size 3 side pads pair - Black Vinyl	XXX
<input type="checkbox"/>	151-6317-05	Size 3 side pads pair - Ebony Grey	XXX
<input type="checkbox"/>	151-6116-05	Size 1 side pads pair - Ebony Grey (White Metal)	XXX
<input type="checkbox"/>	151-6223-05	Size 2 side pads pair - Ebony Grey (White Metal)	XXX
<input type="checkbox"/>	151-6318-05	Size 3 side pads pair - Ebony Grey (White Metal)	XXX

Protraction Pad



<input type="checkbox"/>	151-6218-05VIN	Size 1/2 Protraction pad pair Ebony Grey Vinyl		XXX
<input type="checkbox"/>	151-6218-05	Size 1/2 Protraction pad pair - Ebony Grey		XXX
<input type="checkbox"/>	151-6314-05VIN	Size 3 Protraction Pad pair - Ebony Grey Vinyl		XXX
<input type="checkbox"/>	151-6314-05	Size 3 Protraction Pad pair - Ebony Grey		XXX
<input type="checkbox"/>	151-6219-05	Size 1/2 Protraction pad pair Ebony Grey Vinyl (White Metal)		XXX
<input type="checkbox"/>	151-6315-05	Size 3 Protraction Pad pair - Ebony Grey (White Metal)		XXX

Abduction Pad



<input type="checkbox"/>	151-6224-07VIN	Size 2 Abduction Side Pad - Black Vinyl		XXX
<input type="checkbox"/>	151-6224-05	Size 2 Abduction Side Pad - Ebony Grey		XXX
<input type="checkbox"/>	151-6320-07VIN	Size 3 Abduction Side Pad - Black Vinyl		XXX
<input type="checkbox"/>	151-6320-05	Size 3 Abduction Side Pad - Ebony Grey		XXX
<input type="checkbox"/>	151-6225-05	Size 2 Abduction Side Pad - Ebony Grey (White Metal)		XXX
<input type="checkbox"/>	151-6321-05	Size 3 Abduction Side Pad - Ebony Grey (White Metal)		XXX

Chest Support



<input type="checkbox"/>	137-808-05	Oval Chest Harness - Black		XXX
<input type="checkbox"/>	137-808-05VIN	Oval Chest Harness - Black Vinyl		XXX
<input type="checkbox"/>	133-625	Oval Chest Harness - Size 3 Only		XXX
<input type="checkbox"/>	155-N001-L600	Small Butterfly Harness & attachment straps		XXX
<input type="checkbox"/>	155-N001-L601	Medium Butterfly Harness & attachment straps		XXX
<input type="checkbox"/>	155-N001-L602	Large Butterfly Harness & attachment straps		XXX
<input type="checkbox"/>	155-N001-L603	X-Large Butterfly Harness & attachment straps		XXX
<input type="checkbox"/>	133-625	EDS size 3 Chest Pad		XXX

Foot Support



Sandals



Booties



Ankle Huggers



Ankle Wraps

<input type="checkbox"/>	152-1600	Sandals Size 1		XXX
<input type="checkbox"/>	152-2600	Sandals Size 2		XXX
<input type="checkbox"/>	152-3600	Sandals Size 3		XXX
<input type="checkbox"/>	150-1852	Extra Small Booties		XXX
<input type="checkbox"/>	AH01	Ankle Huggers		XXX
<input type="checkbox"/>	173-2305	Small Ankle Wraps (use with sandals)		XXX
<input type="checkbox"/>	173-3305	Large Ankle Wraps (use with sandals)		XXX
<input type="checkbox"/>	117-787	1 Inch Sandal Raisers (pair) for split footplate		XXX

Pelvic Support



<input type="checkbox"/>	151-1660	4 point pelvic belt - Size 1		XXX
<input type="checkbox"/>	151-2660	4 point pelvic belt - Size 2		XXX
<input type="checkbox"/>	151-3660	4 point pelvic belt - Size 3		XXX

Technical Details			
	Size 1	Size 2	Size 3
Maximum load	30kgs 66lbs	50kg 110lbs	100kgs 220lbs
Weight - seat and hi-low chassis	29.5kg 65lbs	34kg 75lbs	49kg 108lbs
Seat height – hi-low chassis	250-600 9.8-23.6	320-600 12.6-23.6	450-600 17.7-23.6
Seat width	160-280 6.3-11	260-380 10.2-15	340-500 13.4-19.7
Seat depth (standard)	150-280 5.9-11	230-380 9-15	330-480 13-1109
Seat depth (extended)	150-280 5.9-11	300-450 11.8-17.7	330-480 13-1109
Lateral chest distance	150-220 5.9-10.7	170-240 6.7-9.4	240-370 9.4-14.6
Back height	300-400 11.8-15.7	400-520 15.7-20.5	450-630 17.7-24.8
Back angle	30° supine 10° prone	30° supine 10° prone	30° supine 10° prone
Armrest height	150-250 5.9-9.8	200-300 7.9-11.8	180-325 7-12.7
Lower leg length	150-250 5.9-9.8	250-450 9.8-17.7	250-450 9.8-17.7
Lower leg angle	-10° to +50°	-10° to +50°	-10° to +50°
Footplate angle	-15° to +10°	-15° to +10°	-15° to +10°
Tilt-in-space	30° supine 10° prone	30° supine 10° prone	30° supine 10° prone