



Leckey Mygo Seat CE

Name	
Date:	
Email	
Contact Number	
Invoice Address	
Delivery Address	

Quote Order



Item Code	Description	Quantity	Price
Seat Shell			
<input type="checkbox"/>	117-610-G	Mygo S1 Grey Seat Shell	XXX
<input type="checkbox"/>	137-600-G	Mygo S2 Grey Seat Shell	XXX
<input type="checkbox"/>	185-600-G	Mygo S3 Grey Seat Shell	XXX
<input type="checkbox"/>	117-610	Mygo S1 Black Seat Shell	XXX
<input type="checkbox"/>	137-600	Mygo S2 Black Seat Shell	XXX
<input type="checkbox"/>	185-600	Mygo S3 Black Seat Shell	XXX
Base			
<input type="checkbox"/>	194-700	Mygo Seat Interface for Sunrise Bases	XXX
<input type="checkbox"/>	102-2753-G	Hi-Low Mygo Grey Indoor Manual Chassis - With Height Adjustable Handle	XXX
<input type="checkbox"/>	102-2754-G	Hi-Low Mygo Grey Indoor Powered Chassis - With Height Adjustable Handle	XXX
<input type="checkbox"/>	102-2766-G	Hi-Low Mygo Grey Indoor Hydraulic Chassis - With Height Adjustable Handle	XXX
<input type="checkbox"/>	102-2767-G	Hi-Low Mygo Grey Indoor/Outdoor Manual Chassis - With Height Adjustable Handle	XXX
<input type="checkbox"/>	102-2768-G	Hi-Low Mygo Grey Indoor/Outdoor Hydraulic Chassis - With Height Adjustable Handle	XXX
<input type="checkbox"/>	102-2769-G	Hi-Low Mygo Grey Indoor/Outdoor Powered Chassis - With Height Adjustable Handle	XXX
<input type="checkbox"/>	102-2753	Hi-Low Mygo Black Indoor Manual Chassis - With Push Handle	XXX
<input type="checkbox"/>	102-2754	Hi-Low Mygo Black Indoor Powered Chassis - With Push Handle	XXX
<input type="checkbox"/>	102-2766	Hi-Low Mygo Black Indoor Hydraulic Chassis - With Push Handle	XXX
<input type="checkbox"/>	102-2767	Hi-Low Mygo Black Indoor/Outdoor Manual Chassis - With Push Handle	XXX
<input type="checkbox"/>	102-2768	Hi-Low Mygo Black Indoor/Outdoor Hydraulic Chassis - With Push Handle	XXX
<input type="checkbox"/>	102-2769	Hi-Low Mygo Black Indoor/Outdoor Powered Chassis - With Push Handle	XXX

Covers



<input type="checkbox"/>	117-1851-10	MYGO Seat Covers Size 1 - Black Vinyl (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	117-1851-11	MYGO Seat Covers Size 1 - Teal (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	117-1851-12	MYGO Seat Covers Size 1 - Blue (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	117-1851-13	MYGO Seat Covers Size 1 - Red (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	117-1851-09	MYGO Seat Covers Size 1 - Pink Bubbles (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	117-1851-08	MYGO Seat Covers Size 1 - Blue Bubbles (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	117-1851-07	MYGO Seat Covers Size 1 - Orange Bubbles (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	117-1851-06	MYGO Seat Covers Size 1 - Grey Bubbles (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	117-2851-10	MYGO Seat Covers Size 2 - Black Vinyl (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	117-2851-11	MYGO Seat Covers Size 2 - Teal (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	117-2851-12	MYGO Seat Covers Size 2 - Blue (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	117-2851-13	MYGO Seat Covers Size 2 - Red (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	117-2851-09	MYGO Seat Covers Size 2 - Pink Bubbles (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	117-2851-08	MYGO Seat Covers Size 2 - Blue Bubbles (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	117-2851-07	MYGO Seat Covers Size 2 - Orange Bubbles (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	117-2851-06	MYGO Seat Covers Size 2 - Grey Bubbles (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	185-2851-10	MYGO Seat Covers Size 3 - Black Vinyl (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	185-2851-11	MYGO Seat Covers Size 3 - Teal (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	185-2851-12	MYGO Seat Covers Size 3 - Blue (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	185-2851-13	MYGO Seat Covers Size 3 - Red (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	185-2851-09	MYGO Seat Covers Size 3 - Pink Bubbles (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	185-2851-08	MYGO Seat Covers Size 3 - Blue Bubbles (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	185-2851-07	MYGO Seat Covers Size 3 - Orange Bubbles (includes seat, sacral and backrest)	XXX
<input type="checkbox"/>	185-2851-06	MYGO Seat Covers Size 3 - Grey Bubbles (includes seat, sacral and backrest)	XXX

Footplate



Single Footplate



Split Footplate

<input type="checkbox"/>	117-3776	Single footplate - inc padded cover Size 1	XXX
<input type="checkbox"/>	117-4776	Single footplate - inc padded cover Size 2	XXX
<input type="checkbox"/>	117-1776	Split Footplate Size 1	XXX
<input type="checkbox"/>	117-2776	Split Footplate Size 2	XXX

Hip Guides



Size 0



Size 1/2

<input type="checkbox"/>	137-632	MYGO S0 Hip Laterals		XXX
<input type="checkbox"/>	137-632-05VIN	MYGO S0 Hip Laterals - Vinyl		XXX
<input type="checkbox"/>	117-642-05	Mygo S1 Hip Guide Size and covers		XXX
<input type="checkbox"/>	117-642-05VIN	Mygo S1 Hip Guide Size and vinyl covers		XXX
<input type="checkbox"/>	117-643-05	Mygo S2 Hip Guide Size and covers		XXX
<input type="checkbox"/>	117-643-05VIN	Mygo S2 Hip Guide Size and vinyl covers		XXX

Armrest




Standard



Wide

<input type="checkbox"/>	117-1777	MYGO Arm Rest Size 1		XXX
<input type="checkbox"/>	117-2777	MYGO Arm Rest Size 2		XXX
<input type="checkbox"/>	185-604	MYGO Arm Rest Size 3		XXX
<input type="checkbox"/>	117-645	Mygo Size 1/2 Wide armrest (pair)		XXX


Tray



Black Tray



Grab Rails



Padded Tray



Elbow Blockers

<input type="checkbox"/>	117-646	Mygo S1&2 Black Tray		XXX
<input type="checkbox"/>	117-647	Mygo S3 Black Tray		XXX
<input type="checkbox"/>	117-652	Black Padded Tray (Size 1 and 2)		XXX
<input type="checkbox"/>	117-641	Pair Tray elbow blockers inc. cover		XXX
<input type="checkbox"/>	117-769	Grab Rails		XXX

Spacer Pad



<input type="checkbox"/>	117-768-09	MYGO Pink 1" Back Spacer Pad		XXX
<input type="checkbox"/>	117-768-08	MYGO Blue 1" Back Spacer Pad		XXX
<input type="checkbox"/>	117-768-07	MYGO Orange 1" Back Spacer Pad		XXX
<input type="checkbox"/>	117-768-06	MYGO Grey 1" Back Spacer Pad		XXX
<input type="checkbox"/>	117-768-05	Mygo Black Back Spacer Pad		XXX
<input type="checkbox"/>	117-768-05VIN	Mygo Black Vinyl 1" Back Spacer Pad		XXX

Shoulder Support



<input type="checkbox"/>	137-621	MYGO Shoulder Support Hardware		XXX
<input type="checkbox"/>	137-801-09	MYGO Pink Bubbles Shoulder Support Cushion		XXX
<input type="checkbox"/>	137-801-08	MYGO Blue Bubbles Shoulder Support Cushion		XXX
<input type="checkbox"/>	137-801-07	MYGO Orange Bubbles Shoulder Support Cushion		XXX
<input type="checkbox"/>	137-801-06	MYGO Grey Bubbles Shoulder Support Cushion		XXX
<input type="checkbox"/>	137-801-10	MYGO Black Vinyl Shoulder Support Cushion		XXX
<input type="checkbox"/>	137-801-11	MYGO Teal Shoulder Support Cushion		XXX
<input type="checkbox"/>	137-801-12	MYGO New Blue Shoulder Support Cushion		XXX
<input type="checkbox"/>	137-801-13	MYGO Red Shoulder Support Cushion		XXX

Dynamic Backrest



<input type="checkbox"/>	137-662	Dynamic Backrest 150N		XXX
<input type="checkbox"/>	137-661	Dynamic Backrest 100N		XXX

Accessories



<input type="checkbox"/>	117-765-05VIN	MYGO Knee Pads (Vinyl)		XXX
<input type="checkbox"/>	117-847-10	MYGO BLACK VINYL One Piece Cushion Expander		XXX

Chest Support



<input type="checkbox"/>	137-808-05	Oval Chest Harness - Black		XXX
<input type="checkbox"/>	137-808-05VIN	Oval Chest Harness - Black Vinyl		XXX
<input type="checkbox"/>	155-N001-L600	Small Butterfly Harness & attachment straps		XXX
<input type="checkbox"/>	155-N001-L601	Medium Butterfly Harness & attachment straps		XXX
<input type="checkbox"/>	155-N001-L602	Large Butterfly Harness & attachment straps		XXX
<input type="checkbox"/>	155-N001-L603	X-Large Butterfly Harness & attachment straps		XXX

Foot Support



Sandals



Ankle Huggers



Ankle Wraps

<input type="checkbox"/>	152-1600	Sandals Size 1		XXX
<input type="checkbox"/>	152-2600	Sandals Size 2		XXX
<input type="checkbox"/>	152-3600	Sandals Size 3		XXX
<input type="checkbox"/>	AH01	Ankle Huggers		XXX
<input type="checkbox"/>	173-2305	Small Ankle Wraps (use with sandals)		XXX
<input type="checkbox"/>	173-3305	Large Ankle Wraps (use with sandals)		XXX
<input type="checkbox"/>	117-787	1Inch Sandal Raisers (pair) for split footplate		XXX

Headrest



Contoured



Contoured Occipital



Flat



Laterals (for use with flat)

<input type="checkbox"/>	137-699-05	Contoured headrest Assembly and Cushion		XXX
<input type="checkbox"/>	137-699-05VIN	Contoured headrest Assembly and Cushion - Vinyl		XXX
<input type="checkbox"/>	117-699-05	Contoured Occipital Headrest Foam & Cover		XXX
<input type="checkbox"/>	117-699-05VIN	Contoured Occipital Headrest Foam & Cover - Vinyl		XXX
<input type="checkbox"/>	117-851-05	Flat headrest Assembly and Cushion		XXX
<input type="checkbox"/>	117-851-05VIN	Flat headrest Assembly and Cushion - Vinyl		XXX
<input type="checkbox"/>	120-893-05	Flat Headrest Laterals with covers		XXX
<input type="checkbox"/>	120-893-05VIN	Flat Headrest Laterals with covers Vinyl		XXX
<input type="checkbox"/>	117-449	1" Wooden spacer for 8 Inch PRO Bar		XXX

Pelvic Support



2 Point Belt



Pelvic Cradle



Pelvic Harness



Pelvic Spacers

<input type="checkbox"/>	117-1860	2-point belt - Size 1		XXX
<input type="checkbox"/>	117-2860	2-point belt - Size 2		XXX
<input type="checkbox"/>	117-3860	2-point belt - Size 3		XXX
<input type="checkbox"/>	117-476	Pelvic Cradle Size 2 (old size 00)		XXX
<input type="checkbox"/>	117-477	Pelvic Cradle Size 3 (old size 0)		XXX
<input type="checkbox"/>	117-478	Pelvic cradle Size 4 (old size 1)		XXX
<input type="checkbox"/>	117-479	Pelvic cradle Size 5 (old size 2)		XXX
<input type="checkbox"/>	117-480	Pelvic cradle Size 6 (old size 3)		XXX
<input type="checkbox"/>	133-737	Pelvic Cradle Gluteal Ramp (Size 0-3)		XXX
<input type="checkbox"/>	117-802	MYGO Pelvic Harness & Brackets - Medium		XXX
<input type="checkbox"/>	117-801	MYGO Pelvic Harness & Brackets - Small		XXX
<input type="checkbox"/>	120-789	MYGO/SQSE Pelvic Harness spacer pads (pair)		XXX

Laterals



Rigid Laterals



Flipaway Laterals



Laterals with Protraction Pads



Complex Laterals

<input type="checkbox"/>	117-803	Rigid Laterals (with cover)		XXX
<input type="checkbox"/>	117-803-05VIN	Vinyl Rigid Laterals (with cover)		XXX
<input type="checkbox"/>	117-831	Spacer Laterals		XXX
<input type="checkbox"/>	117-831-05VIN	Black Vinyl Spacer Laterals		XXX
<input type="checkbox"/>	117-838	Narrow Laterals (1" narrower chest width)		XXX
<input type="checkbox"/>	117-838-05VIN	Black Narrow Laterals (1" narrower chest width)		XXX
<input type="checkbox"/>	117-837	Flipaway Lateral Supports		XXX
<input type="checkbox"/>	117-837-05VIN	Black Vinyl Flipaway Lateral Supports		XXX
<input type="checkbox"/>	117-644-05	Soft Flipaway laterals with protraction pad - pair with covers		XXX
<input type="checkbox"/>	117-644-05VIN	Black Vinyl Flipaway laterals with protraction pad - pair with covers		XXX
<input type="checkbox"/>	137-634	Flipaway Laterals Extended Height		XXX
<input type="checkbox"/>	137-634-05VIN	Black Vinyl Flipaway Laterals Extended Height		XXX
<input type="checkbox"/>	137-633	Rigid Laterals Extended Height 50 mm - pair		XXX
<input type="checkbox"/>	137-633-05VIN	Black Vinyl Rigid Laterals Extended Height 50 mm - pair		XXX
<input type="checkbox"/>	137-193	Rigid Laterals Extended Height 100 mm - pair		XXX
<input type="checkbox"/>	137-193-05VIN	Black Vinyl Rigid Laterals Extended Height 100 mm - pair		XXX
<input type="checkbox"/>	137-638	Right Hand Horizontal Complex Lateral		XXX
<input type="checkbox"/>	137-637	Left Hand Horizontal Complex Lateral		XXX

	Size 1	Size 2	Size 3
Age (approx)	3-10	8-14	10-16
User Weight	Min 18kg [40lbs] Max 50kg [110lbs]	Min 18kg [40lbs] Max 60kg [132lbs]	Min 18kg [40lbs] Max 70kg [154lbs]
User Height	Min 105cm [41"] Max 150cm [59"]	Min 127cm [50"] Max 168cm [66"]	Min 150cm [59"] Max 175cm [69"]
Seat Width	Min 200mm [8"] Max 325mm [13"]	Min 220mm [8.7"] Max 345mm [13.6"]	Min 250mm [9.8"] Max 400mm [15.7"]
Seat Depth	Min 270mm [10.6"] Max 420mm [16.5"]	Min 350mm [13.8"] Max 470mm [18.5"]	Min 350mm [13.8"] Max 470mm [18.5"]
Chest Width	Min 170mm [6.7"] Max 270mm [10.6"]	Min 170mm [6.7"] Max 270mm [10.6"]	Min 200mm [7.8"] Max 300mm [11.8"]
Backrest Height	Min 360mm [14.2"] Max 470mm [18.5"]	Min 460mm [18.1"] Max 570mm [22.4"]	Min 460mm [18.1"] Max 570mm [22.4"]
Backrest Angle	Prone 10° - Recline 25°	Prone 10° - Recline 25°	Prone 10° - Recline 25°
Lower Leg Length	Min 215mm [8.5"] Max 350mm [13.8"]	Min 315mm [12.4"] Max 470mm [18.5"]	Min 315mm [12.4"] Max 470mm [18.5"]
Armrest Height	Min 160mm [6.3"] Max 210mm [8.3"]	Min 210mm [8.3"] Max 260mm [10.2"]	Min 210mm [8.3"] Max 260mm [10.2"]
Tray Size	550 x 480mm [21.6" x 18.9"]	550 x 480mm [21.6" x 18.9"]	550 x 480mm [21.6" x 18.9"]
Seat Height (Footpedal Chassis)	340mm-655mm [13.4"-25.8"]	340mm-655mm [13.4"-25.8"]	340mm-655mm [13.4"-25.8"]
Seat Height (Powered Chassis)	370mm-675mm [14.6"-26.6"]	370mm-675mm [14.6"-26.6"]	370mm-675mm [14.6"-26.6"]