



# Leckey PAL CE

Name	
Date:	
Email	
Contact Number	
Invoice Address	
Delivery Address	

Quote     Order

Purchase Order	
----------------	--



Item Code	Description	Quantity	Price
Seat			
<input type="checkbox"/>	135-1600-03	PAL Size 1 Seat (including lap strap) - Blue	XXX
<input type="checkbox"/>	135-1600-01	PAL Size 1 Seat (including lap strap) - Green	XXX
<input type="checkbox"/>	135-2600-03	PAL Size 2 Seat (including lap strap) - Blue	XXX
<input type="checkbox"/>	135-2600-01	PAL Size 2 Seat (including lap strap) - Green	XXX
<input type="checkbox"/>	135-3600-03	PAL Size 3 Seat (including lap strap) - Blue	XXX
<input type="checkbox"/>	135-3600-01	PAL Size 3 Seat (including lap strap) - Green	XXX
<input type="checkbox"/>	135-4600-03	PAL Size 4 Seat (including lap strap) - Blue	XXX
<input type="checkbox"/>	135-4600-01	PAL Size 4 Seat (including lap strap) - Green	XXX
Tray			
<input type="checkbox"/>	135-1620	PAL Activity tray size 1	XXX
<input type="checkbox"/>	135-1656	PAL Tray Pad size 1	XXX
<input type="checkbox"/>	135-2620	PAL Activity tray size 2	XXX
<input type="checkbox"/>	135-2656	PAL Tray Pad size 2	XXX
<input type="checkbox"/>	135-3620	PAL Activity tray size 3	XXX
<input type="checkbox"/>	135-3656	PAL Tray Pad size 3	XXX
<input type="checkbox"/>	135-4620	PAL Activity tray size 4	XXX
<input type="checkbox"/>	135-4656	PAL Tray Pad size 4	XXX

### Side Pads



<input type="checkbox"/>	135-1635-03	PAL 25mm Side Pads size 1 (pair) - Blue		XXX
<input type="checkbox"/>	135-1635-01	PAL 25mm Side Pads size 1 (pair) - Green		XXX
<input type="checkbox"/>	135-1654-03	PAL 50mm Side Pads size 1 (pair) - Blue		XXX
<input type="checkbox"/>	135-1654-01	PAL 50mm Side Pads size 1 (pair) - Green		XXX
<input type="checkbox"/>	135-2635-03	PAL 25mm Side Pads size 2 (pair) - Blue		XXX
<input type="checkbox"/>	135-2635-01	PAL 25mm Side Pads size 2 (pair) - Green		XXX
<input type="checkbox"/>	135-2654-03	PAL 50mm Side Pads size 2 (pair) - Blue		XXX
<input type="checkbox"/>	135-2654-01	PAL 50mm Side Pads size 2 (pair) - Green		XXX
<input type="checkbox"/>	135-3635-03	PAL 25mm Side Pads size 3-4 (pair) - Blue		XXX
<input type="checkbox"/>	135-3635-01	PAL 25mm Side Pads size 3-4 (pair) - Green		XXX
<input type="checkbox"/>	135-3654-03	PAL 50mm Side Pads size 3-4 (pair) - Blue		XXX
<input type="checkbox"/>	135-3654-01	PAL 50mm Side Pads size 3-4 (pair) - Green		XXX

### Ramped Cushion



<input type="checkbox"/>	135-1651-03	PAL Ramped Seat Cushion size 1 - Blue		XXX
<input type="checkbox"/>	135-1651-01	PAL Ramped Seat Cushion size 1 - Green		XXX
<input type="checkbox"/>	135-2651-03	PAL Ramped Seat Cushion size 2 - Blue		XXX
<input type="checkbox"/>	135-2651-01	PAL Ramped Seat Cushion size 2 - Green		XXX
<input type="checkbox"/>	135-3651-03	PAL Ramped Seat Cushion size 3 - Blue		XXX
<input type="checkbox"/>	135-3651-01	PAL Ramped Seat Cushion size 3 - Green		XXX
<input type="checkbox"/>	135-4651-03	PAL Ramped Seat Cushion size 4 - Blue		XXX
<input type="checkbox"/>	135-4651-01	PAL Ramped Seat Cushion size 4 - Green		XXX

### Pommel



<input type="checkbox"/>	135-1602-03	PAL Pommel (size 1) - Blue		XXX
<input type="checkbox"/>	135-1602-01	PAL Pommel (size 1) - Green		XXX
<input type="checkbox"/>	135-3602-03	PAL Pommel (size 2 3 4) - Blue		XXX
<input type="checkbox"/>	135-3602-01	PAL Pommel (size 2 3 4) - Green		XXX

### Footrest



<input type="checkbox"/>	135-1603	PAL Footrest size 1		XXX
<input type="checkbox"/>	135-2603	PAL Footrest size 2		XXX
<input type="checkbox"/>	135-3603	PAL Footrest size 3		XXX
<input type="checkbox"/>	135-4603	PAL Footrest size 4		XXX

### Base



Rockers



Stabilisers



Skis with castors

<input type="checkbox"/>	135-1610	PAL Ski with castors size 1 (pair)		XXX
<input type="checkbox"/>	135-1604	PAL Rocker / stabilisers size 1 (pair)		XXX
<input type="checkbox"/>	135-2610	PAL Ski with castors size 2 (pair)		XXX
<input type="checkbox"/>	135-2604	PAL Rocker / stabilisers size 2 (pair)		XXX
<input type="checkbox"/>	135-3610	PAL Ski with castors size 3 (pair)		XXX
<input type="checkbox"/>	135-3604	PAL Rocker / stabilisers size 3 (pair)		XXX
<input type="checkbox"/>	135-4610	PAL Ski with castors size 4 (pair)		XXX
<input type="checkbox"/>	135-4604	PAL Rocker / stabilisers size 4 (pair)		XXX

	Size 1	Size 2	Size 3	Size 4
Age	1-3	2-5	4-8	7-12
User Height	Min 700mm [27.6"] Max 950mm [37.4"]	Min 900mm [35.4"] Max 1050mm [41.3"]	Min 1000mm [39.4"] Max 1300mm [51.2"]	Min 1200mm [47"] Max 1600mm [63"]
Max. User Weight	30kg	35kg	40kg	50kg
Seat Depth	180mm-240mm [7"-9.4"]	220mm-290mm [8.6"-11.4"]	250mm-330mm [9.8"-12.9"]	300mm-400mm [11.8"-15.7"]
Seat Width	150mm-250mm [5.9"-9.8"]	170mm-270mm [6.6"-10.6"]	230mm-330mm [9"-12.9"]	280mm-380mm [11"-14.9"]
Seat Height	160mm-240mm [6.2"-9.4"]	185mm-285mm [7.2"-11.2"]	250mm-350mm [9.8"-13.7"]	300mm-400mm [11.8"-15.7"]