






Leckey Squiggles + CE

Name	
Date:	
Email	
Contact Number	
Invoice Address	
Delivery Address	



Item Code	Description	Quantity	Price
Frame and Chassis			
			
<input type="checkbox"/>	184-764	Squiggle + Pivot Chassis (gas operated - can be adjusted with child in frame)	XXX
<input type="checkbox"/>	184-600	Squiggle + Stander Support Shell (central support includes hip & chest support with sternum prompt chest & hip lateral supports kneecups footplate and tray attachment)	XXX
Cover Pack Includes 2 x support harnesses, 2 x support cushions, 2 x lateral covers, 1 x pair of kneepads			
			
<input type="checkbox"/>	184-1800-09VIN	Cover Pack (Size 1 Knees)- Red Vinyl	XXX
<input type="checkbox"/>	184-1800-09	Cover Pack (Size 1 Knees) - Red	XXX
<input type="checkbox"/>	184-1800-04	Cover Pack (Size 1 Knees)- Purple	XXX
<input type="checkbox"/>	184-1800-02	Cover Pack (Size 1 Knees)- Blue	XXX
<input type="checkbox"/>	184-2800-09VIN	Cover Pack (Size 2 Knees) - Red Vinyl	XXX
<input type="checkbox"/>	184-2800-09	Cover Pack (Size 2 Knees)- Red	XXX
<input type="checkbox"/>	184-2800-04	Cover Pack (Size 2 Knees)- Purple	XXX
<input type="checkbox"/>	184-2800-02	Cover Pack (Size 2 Knees)- Blue	XXX
Tray			
			
<input type="checkbox"/>	118-619	Squiggles Stander Activity Tray	XXX
<input type="checkbox"/>	188-1854-05	Squiggles Black Padded Tray Cover	XXX

Foot Supports				
				
	Sandals	Booties	Ankle Huggers	Ankle Wraps
<input type="checkbox"/>	152-1600	Sandals Size 1		XXX
<input type="checkbox"/>	AH01	Ankle Huggers		XXX
<input type="checkbox"/>	150-1852	Extra Small Booties		XXX
<input type="checkbox"/>	173-2305	Small Ankle Wraps (use with sandals)		XXX
<input type="checkbox"/>	117-787	1Inch Footplate Raisers for Split Footplate		XXX
Head Supports				
				
	Contoured	Contoured Occipital	Flat	Laterals (for use with flat)
<input type="checkbox"/>	137-699-05	Contoured headrest Assembly and Cushion		XXX
<input type="checkbox"/>	137-699-05VIN	Contoured headrest Assembly and Cushion - Vinyl		XXX
<input type="checkbox"/>	117-699-05	Contoured Occipital Headrest		XXX
<input type="checkbox"/>	117-699-05VIN	Contoured Occipital Headrest - Vinyl		XXX
<input type="checkbox"/>	117-851-05	Flat Headrest Assembly and Cushion		XXX
<input type="checkbox"/>	117-851-05VIN	Flat Headrest Assembly and Cushion - Vinyl		XXX
<input type="checkbox"/>	120-893-05	Flat Headrest Laterals with covers		XXX
<input type="checkbox"/>	120-893-05VIN	Flat Headrest Laterals with covers Vinyl		XXX
Targetted Training Conversion Kit (includes dynamic base and wraparound supports)				
				
<input type="checkbox"/>	188-680-02	Squiggle + to Squiggle TT conversion kit blue		XXX
<input type="checkbox"/>	188-680-04	Squiggle + to Squiggle TT conversion kit purple		XXX
<input type="checkbox"/>	188-680-09	Squiggle + to Squiggle TT conversion kit red		XXX
<input type="checkbox"/>	188-680-09VIN	Squiggle + to Squiggle TT conversion kit red vinyl		XXX

Age (approx)	1-5
Max User Weight	22kg/48.4lbs
User Height	Min 730mm / 28.5 inches Max 1110mm / 43.5 inches
Chest Support	Angle adjustment +/- 30° Depth adjustment 50mm / 2 inches
Chest Support Height (footplate to armpit)	Min 460mm / 18.1 inches Max 820mm / 32 inches
Chest Width (distance between laterals)	Min 160mm / 6.3 inches Max 230mm / 9 inches
Hip Support (height from footplate to mid buttock)	Min 300mm / 11.8 inches Max 530mm / 20.8 inches
Distance between midline of knees	Min 140mm / 5.5 inches Max 215mm / 8.4 inches
Footplate Angle	Plantarflexion 10° / Dorsiflexion 10°