







Leckey Squiggles Seat CE

Name	
Date:	
Email	
Contact Number	
Invoice Address	
Delivery Address	

Quote Order

Purchase Order	
----------------	--



Item Code	Description	Quantity	Price
Seat			
			
<input type="checkbox"/>	120-610-MIN	Squiggles Seat Squiggles Shell-Configured in Minimum Seat Depth Set-up	XXX
Interface			
			
<input type="checkbox"/>	194-600	Squiggles Seat Interface for Sunrise bases	XXX
High Low Base			
 			
<input type="checkbox"/>	120-800	Squiggles Seat Hi low chassis - foot pedal	XXX
<input type="checkbox"/>	122-710	Mygo/Squiggles Standard Black Push handle	XXX
<input type="checkbox"/>	122-2710	Mygo/Squiggles Height Adjustable Push Handle (Grey Metal)	XXX

Mygo Base



Indoor

Indoor/Outdoor

<input type="checkbox"/>	120-818	Interface to Attach Squiggles Seat to Mygo Chassis		XXX
<input type="checkbox"/>	102-2753-G	Hi-Low Mygo Grey Indoor Manual Chassis - With Height Adjustable Handle		XXX
<input type="checkbox"/>	102-2754-G	Hi-Low Mygo Grey Indoor Powered Chassis - With Height Adjustable Handle		XXX
<input type="checkbox"/>	102-2766-G	Hi-Low Mygo Grey Indoor Hydraulic Chassis - With Height Adjustable Handle		XXX
<input type="checkbox"/>	102-2767-G	Hi-Low Mygo Grey Indoor/Outdoor Manual Chassis - With Height Adjustable Handle		XXX
<input type="checkbox"/>	102-2768-G	Hi-Low Mygo Grey Indoor/Outdoor Hydraulic Chassis - With Height Adjustable Handle		XXX
<input type="checkbox"/>	102-2769-G	Hi-Low Mygo Grey Indoor/Outdoor Powered Chassis - With Height Adjustable Handle		XXX
<input type="checkbox"/>	102-2753	Hi-Low Mygo Black Indoor Manual Chassis - With Push Handle		XXX
<input type="checkbox"/>	102-2754	Hi-Low Mygo Black Indoor Powered Chassis - With Push Handle		XXX
<input type="checkbox"/>	102-2766	Hi-Low Mygo Black Indoor Hydraulic Chassis - With Push Handle		XXX
<input type="checkbox"/>	102-2767	Hi-Low Mygo Black Indoor/Outdoor Manual Chassis - With Push Handle		XXX
<input type="checkbox"/>	102-2768	Hi-Low Mygo Black Indoor/Outdoor Hydraulic Chassis - With Push Handle		XXX
<input type="checkbox"/>	102-2769	Hi-Low Mygo Black Indoor/Outdoor Powered Chassis - With Push Handle		XXX

Tray



Tray




Padded Tray Cover



Grab Rail

<input type="checkbox"/>	120-672	Squiggles Seat tray		XXX
<input type="checkbox"/>	120-688	Squiggles Seat tray 50 mm extended height		XXX
<input type="checkbox"/>	117-769	MYGO/Squiggles Grab Rails		XXX
<input type="checkbox"/>	188-1854-05	Squiggles Black Padded Tray Cover		XXX

Covers



Blue

Green

Pink

<input type="checkbox"/>	120-710-04	Squiggles Seat covers - pink (Includes femorals seat and backrest)		XXX
<input type="checkbox"/>	120-710-03	Squiggles Seat covers - blue (Includes femorals seat and backrest)		XXX
<input type="checkbox"/>	120-710-01	Squiggles Seat covers - green (Includes femorals seat and backrest)		XXX

Shoulder Support



<input type="checkbox"/>	120-702	Squiggles Seat Shoulder Support Assembly		XXX
<input type="checkbox"/>	120-828-04	Squiggles Seat Pink Shoulder Cushion		XXX
<input type="checkbox"/>	120-828-03	Squiggles Seat Blue Shoulder Cushion		XXX
<input type="checkbox"/>	120-828-01	Squiggles Seat Green Shoulder Cushion		XXX

Leg Guide				
				
<input type="checkbox"/>	120-856-04	Squiggles Seat Leg Guide Pad Assembly - pink		XXX
<input type="checkbox"/>	120-856-03	Squiggles Seat Leg Guide Pad Assembly - blue		XXX
<input type="checkbox"/>	120-856-01	Squiggles Seat Leg Guide Pad Assembly - green		XXX
Footrest and Pommel				
 				
<input type="checkbox"/>	120-641	Squiggles Seat pommel		XXX
<input type="checkbox"/>	120-639	Squiggles Seat Footrest		XXX
Back Spacer				
				
<input type="checkbox"/>	120-779-04	Squiggles Seat 1" back spacer pad & flipaway lateral support - Pink		XXX
<input type="checkbox"/>	120-779-03	Squiggles Seat 1" back spacer pad & flipaway lateral support - Blue		XXX
<input type="checkbox"/>	120-779-01	Squiggles Seat 1" back spacer pad & flipaway lateral support - Green		XXX
<input type="checkbox"/>	120-777-04	Squiggles Seat 1" back spacer pad & rigid lateral support - Pink		XXX
<input type="checkbox"/>	120-777-03	Squiggles Seat 1" back spacer pad & rigid lateral support - Blue		XXX
<input type="checkbox"/>	120-777-01	Squiggles Seat 1" back spacer pad & rigid lateral support - Green		XXX
Chest Support				
 				
<input type="checkbox"/>	137-808-05	Oval Chest Harness - Black		XXX
<input type="checkbox"/>	137-808-05VIN	Oval Chest Harness - Black Vinyl		XXX
<input type="checkbox"/>	155-N001-L600	Small Butterfly Harness & attachment straps		XXX
<input type="checkbox"/>	155-N001-L601	Medium Butterfly Harness & attachment straps		XXX
Foot Support				
   				
		Sandals		
		Booties		
		Ankle Huggers		
		Ankle Wraps		
<input type="checkbox"/>	152-1600	Sandals Size 1		XXX
<input type="checkbox"/>	150-1852	Extra Small Booties		XXX
<input type="checkbox"/>	AH01	Ankle Huggers		XXX
<input type="checkbox"/>	173-2305	Small Ankle Wraps (use with sandals)		XXX
<input type="checkbox"/>	120-816	1Inch Footplate Raiser for Single Footplate		XXX
Head Support				



Contoured



Contoured Occipital



Flat



Laterals (for use with flat)

<input type="checkbox"/>	137-699-05	Contoured headrest Assembly and Cushion		XXX
<input type="checkbox"/>	137-699-05VIN	Contoured headrest Assembly and Cushion - Vinyl		XXX
<input type="checkbox"/>	117-699-05	Contoured Occipital Headrest		XXX
<input type="checkbox"/>	117-699-05VIN	Contoured Occipital Headrest - Vinyl		XXX
<input type="checkbox"/>	117-851-05	Flat Headrest Assembly and Cushion		XXX
<input type="checkbox"/>	117-851-05VIN	Flat Headrest Assembly and Cushion - Vinyl		XXX
<input type="checkbox"/>	120-893-05	Flat Headrest Laterals with covers		XXX
<input type="checkbox"/>	120-893-05VIN	Flat Headrest Laterals with covers Vinyl		XXX

Pelvic Support



Pelvic Harness



Pelvic Spacer

<input type="checkbox"/>	120-844-03	Squiggles pelvic harness small - blue		XXX
<input type="checkbox"/>	120-844-01	Squiggles pelvic harness small - green		XXX
<input type="checkbox"/>	120-844-04	Squiggles pelvic harness small - pink		XXX
<input type="checkbox"/>	117-848-03	Squiggles pelvic harness medium - blue		XXX
<input type="checkbox"/>	117-848-01	Squiggles pelvic harness medium - green		XXX
<input type="checkbox"/>	117-848-04	Squiggles pelvic harness medium - pink		XXX
<input type="checkbox"/>	120-789	MYGO/Squiggles Pelvic Harness spacer pads (pair)		XXX

Laterals



Rigid Laterals



Flipaway Laterals



Flipaway Laterals with Protraction Pads

<input type="checkbox"/>	120-656	Squiggles Seat rigid lateral support - Hardware		XXX
<input type="checkbox"/>	120-855	Squiggles Seat flipaway lateral support Hardware		XXX
<input type="checkbox"/>	120-838-03	Squiggles Seat lateral support covers - blue		XXX
<input type="checkbox"/>	120-838-01	Squiggles Seat lateral support covers - green		XXX
<input type="checkbox"/>	120-838-04	Squiggles Seat lateral support covers - pink		XXX
<input type="checkbox"/>	117-644-05	Flipaway Laterals with Protraction Pads		XXX

Product Sizing			
		mm	Inches
Age (Approx)	1-5 Years		
Stature	Min	750	29.5
	Max	1100	43.3
Max User Weight (kg/Lbs)		22Kg	48.4Lbs
Seat Depth (Backrest to back of knee)	Min	190	7.5
	Max	290	11.4
Seat width	Min	160	6.3
	Max	260	10.2
Hip Width (Distance between hip guides)	Min	160	6.3
	Max	220	8.7
Backrest Height (Sitting shoulder height)	Min	320	12.6
	Max	380	15
Backrest Height (with shoulder support)	Min	420	16.5
	Max	480	18.9
Chest Width (Distance between laterals)	Min	160	6.3
	Max	90	3.5
Knee Width	Min	190	7.5
	Max	122.5	4.8
Distance between pommel & side of seat Seat to Footplate (Sole of foot to back of knee)	Min	270	10
	Max	650	25
Seat to Floor (Top of Seat to Floor (Hi-Low Chassis))	Prone	10°	10°
	Recline	25°	25°
Footplate Angle	Plantarflexion	10°	10°
	Dorsiflexion	10°	10°
Tray Size		460/560	18.1/22