



Leckey Squiggles Stander

CE

Name	
Date:	
Email	
Contact Number	
Invoice Address	
Delivery Address	

Quote Order

Purchase Order	
----------------	--



	Item Code	Description	Quantity	Price
Frame and Chassis				
				
<input type="checkbox"/>	118-764	Squiggles Stander Pivot Chassis (gas operated - can be adjusted with child in frame)		XXX
<input type="checkbox"/>	118-600	Squiggles Stander Stander Support Shell (central support includes hip & chest support with sternum prompt chest & hip lateral supports kneecups footplate and tray attachment)		XXX
Covers				
				
<input type="checkbox"/>	118-763	Squiggles Stander Stander Covers - Pink		XXX
<input type="checkbox"/>	118-762-03VIN	Squiggles Stander Stander Covers - Blue Vinyl		XXX
<input type="checkbox"/>	118-762	Squiggles Stander Stander Covers - Blue		XXX
<input type="checkbox"/>	118-760	Squiggles Stander Stander Covers - Green		XXX
Tray				
				
				
<input type="checkbox"/>	118-619	Squiggles Stander Activity Tray		XXX
<input type="checkbox"/>	188-1854-05	Squiggles Black Padded Tray Cover		XXX
Foot Supports				
				
<input type="checkbox"/>	152-1600	Sandals Size 1		XXX
<input type="checkbox"/>	150-1852	Extra Small Booties		XXX
<input type="checkbox"/>	AH01	Ankle Huggers		XXX
<input type="checkbox"/>	173-2305	Small Ankle Wraps (use with sandals)		XXX
<input type="checkbox"/>	120-816	1Inch Footplate Raiser for Single Footplate		XXX

Head Supports



Contoured



Contoured Occipital



Flat



Laterals (for use with flat)

<input type="checkbox"/>	137-699-05	Contoured headrest Assembly and Cushion		XXX
<input type="checkbox"/>	137-699-05VIN	Contoured headrest Assembly and Cushion - Vinyl		XXX
<input type="checkbox"/>	117-699-05	Contoured Occipital Headrest		XXX
<input type="checkbox"/>	117-699-05VIN	Contoured Occipital Headrest - Vinyl		XXX
<input type="checkbox"/>	117-851-05	Flat Headrest Assembly and Cushion		XXX
<input type="checkbox"/>	117-851-05VIN	Flat Headrest Assembly and Cushion - Vinyl		XXX
<input type="checkbox"/>	120-893-05	Flat Headrest Laterals with covers		XXX
<input type="checkbox"/>	120-893-05VIN	Flat Headrest Laterals with covers Vinyl		XXX

	One Size Only
Age	1-5
User Height	730mm-1100mm [28.7"-43.3"]
Max. User Weight	22kg (48lbs)
Chest Support Height (footplate to armpit)	550mm-820mm [21.6"-32.2"]
Chest Width (distance between laterals)	160mm-230mm [6.2"-9"]
Chest Support Adjustment	+/- 30° 50mm [1.9"] - Depth
Hip Support (height from footplate to midbuttocks)	300mm-530mm [11.8"-20.8"]
Hip Width (distance between hip guides)	160mm-230mm [6.2"-9"]
Distance between midline of knees	140mm-215mm [5.5"-8.4"]
Footplate Angle	Plantarflexion 10° Dorsiflexion 10°
Tray Height (footplate to top of tray)	555mm-825mm [21.8"-32.4"] - Prone 480mm-825mm [18.8"-32.4"] - Supine
Tray Angle Adjustment	30° - Prone 40° - Supine