





Leckey Therapy Bench CE

Name	
Date:	
Email	
Contact Number	
Invoice Address	
Delivery Address	

Quote Order

Purchase Order	
----------------	--



Item Code	Description	Quantity	Price
Bench			
			
<input type="checkbox"/>	188-660	Leckey Therapy Bench	XXX
Activity Cradles			
			
<input type="checkbox"/>	188-1811	Leckey Therapy Bench size 1 activity cradle	XXX
<input type="checkbox"/>	188-1812	Leckey Therapy Bench size 2 activity cradle	XXX
<input type="checkbox"/>	188-1813	Leckey Therapy Bench size 3 activity cradle	XXX
<input type="checkbox"/>	188-1814	Leckey Therapy Bench size 4 activity cradle	XXX
<input type="checkbox"/>	188-1815	Leckey Therapy Bench size 5 activity cradle	XXX

Foot Support



Sandals



Booties



Ankle Huggers



Ankle Wraps

<input type="checkbox"/>	152-1600	Sandals Size 1		XXX
<input type="checkbox"/>	152-2600	Sandals Size 2		XXX
<input type="checkbox"/>	152-3600	Sandals Size 3		XXX
<input type="checkbox"/>	150-1852	Extra Small Booties		XXX
<input type="checkbox"/>	AH01	Ankle Huggers		XXX
<input type="checkbox"/>	173-2305	Small Ankle Wraps (use with sandals)		XXX
<input type="checkbox"/>	173-3305	Large Ankle Wraps (use with sandals)		XXX

Therapy Bench	
Max User Weight	80kg
Age	1-14
Bench Footprint	620mm x 740mm
Seat Base Height (from ground)	Min 330mm Max 460mm
Seat Base Width	620mm
Seat Base Depth	360mm
Footplate Height (from ground)	Min 20mm Max 110mm
Seat Base to Footplate Height	Min 220mm Max 440mm

Activity Cradle Size Guide		
Cradle Size	Hip Circumference (mm)	ASIS - ASIS (mm)
Size 1	Min 400mm Max 500mm	Under 120
Size 2	Min 440mm Max 575mm	Min 130mm Max 160mm
Size 3	Min 520mm Max 650mm	Min 160mm Max 180mm
Size 4	Min 600mm Max 740mm	Min 190mm Max 200mm
Size 5	Min 660mm Max 855mm	Over 200