



Technical / News Bulletin
UK

BULLETIN NO.:

TB UK209

DATE ISSUED

22/03/2023

BULLETIN DETAILS

Product Group:	JCM	
Model Affected:	Triton	
Effective Date: See Below	Serial No. Effective From:	23/02/2023
Bulletin Subject:	JCM Triton – Pommel Update	

We are committed at Leckey to continually develop products based on customer feedback. Therefore, from listening to customer feedback we are delighted to introduce a new item, the Triton 0 pommel with size 1 bracket.

Previously, customers who wanted to fit a Triton size 0 pommel on to the Triton size 1 could do this however there were significant limitations due to the limited range of adjustment on the bracket.

Therefore, due to customer feedback we have introduced the Size 0 pommel with Size 1 bracketry to allow a greater range of adjustment to ensure better positioning.

The new item can be ordered by using the below part number:

Model	Part Code	Description
Triton	97002828	Triton 0 Pommel with size 1 bracket

For any queries on any of these updates please contact your local Product Specialist.

	Telephone	E-mail
Customer Services	+44 (0)28 9260 0750	hello@leckey.com