



# Leckey Pronestander

<b>Name</b>	
<b>Date:</b>	
<b>Email</b>	
<b>Contact Number</b>	
<b>Invoice Address</b>	
<b>Delivery Address</b>	

Quote     Order

Purchase Order	
----------------	--



	Item Code	Description	Quantity	Price
<b>Frame</b>				
<input type="checkbox"/>	LPS/1-MW	Size 1 Pronestander - Meadow		XXX
<input type="checkbox"/>	LPS/2-MB	Size 2 Pronestander - Mondeo Blue		XXX
<input type="checkbox"/>	LPS/3-MB	Size 3 Pronestander - Mondeo Blue		XXX
<b>Pelvic Band</b>				
<input type="checkbox"/>	LPS/1/B-MW	Pronestander Pelvic Band - size 1 - Meadow		XXX
<input type="checkbox"/>	LPS/2/B-MB	Pronestander Pelvic Band - size 2 - Mondeo Blue		XXX
<input type="checkbox"/>	LPS/3/B-MB	Pronestander Pelvic Band - size 3 - Mondeo Blue		XXX
<input type="checkbox"/>	LFS/1/B-MW	Freestander Pelvic Band - size 1 - Meadow		XXX
<input type="checkbox"/>	LFS/2/B-MB	Freestander Pelvic Band - size 2 - Mondeo Blue		XXX
<input type="checkbox"/>	LFS/3/B-MB	Freestander Pelvic Band - size 3 - Mondeo Blue		XXX
<b>Soft Knees</b>				
<input type="checkbox"/>	LPS/1/K-MW	Freestander/Pronestander Upholstered Kneecups - size 1 - Meadow		XXX
<input type="checkbox"/>	LPS/2/K-MB	Freestander/Pronestander Upholstered Kneecups - size 2 - Mondeo Blue		XXX
<input type="checkbox"/>	LPS/3/K-MB	Freestander/Pronestander Upholstered Kneecups - size 3 - Mondeo Blue		XXX
<input type="checkbox"/>	LPS/4/K-MB	Freestander/Pronestander Upholstered Kneecups - size 4 - Mondeo Blue		XXX
<b>Foot Supports</b>				
<input type="checkbox"/>	152-1600	Sandals Size 1		XXX
<input type="checkbox"/>	152-2600	Sandals Size 2		XXX
<input type="checkbox"/>	152-3600	Sandals Size 3		XXX
<input type="checkbox"/>	LHZ/1/100	Size 1 open-back sandals		XXX
<input type="checkbox"/>	LHZ/2/100	Size 2 open-back sandals		XXX
<input type="checkbox"/>	LHZ/3/100	Size 3 open-back sandals		XXX
<input type="checkbox"/>	AH01	Ankle Huggers		XXX
<input type="checkbox"/>	173-2305	Small Ankle Wraps (use with sandals)		XXX
<input type="checkbox"/>	173-3305	Large Ankle Wraps (use with sandals)		XXX
<input type="checkbox"/>	117-787	1Inch Footplate Raiser for Split Footplate		XXX

Technical data	Size 1	Size 2	Size 3
Age (years)	1.5 – 5	4 – 10	9 – 18
Body height (cm/inches)	Min. 80/31.5 Max. 115/45.3	Min. 110/43.3 Max. 145/57	Min. 140/55.1 Max. 185/72.8
User weight (kg/lbs)	60/132	80/176	90/198
Chest pad height (cm/inches)	Min. 47/18.5 Max. 68/26.8	Min. 68/26.8 Max. 98/38.4	Min. 97/38.4 Max. 114/53.1
Width between chest supports (cm/inches)	Min. 15/6 Max. 25/10	Min. 20/8 Max. 31/12	Min. 25/10 Max. 36/14