



Horizon Stander De-Rotational Pelvic Bands

Product Code	Description	Inner Fabric (mm)	Fabric Around Frame (mm)	Hip Circumference (mm)
LHZ/1/181-GRB	Supine Pelvic Band Size 1	990 x 100	1530 x 100	500 - 800
LHZ/2/181-GRB	Supine Pelvic Band Size 2	1230 x 100	1680 x 100	750 - 950
LHZ/3/181-GRB	Supine Pelvic Band Size 3	1440 x 100	1680 x 100	900 - 1150





Mygo Stander De-Rotational Pelvic Bands

Product Code	Description	Inner Fabric (mm)	Fabric Around Frame (mm)	Hip Circumference (mm)
143-1783-XX	Supine Pelvic Band Size 1	730 x 203	1276 x 150	370 - 590
143-2783-XX	Supine Pelvic Band Size 2	902 x 241	1384 x 191	460 - 730
143-3783-XX	Supine Pelvic Band Size 3	1086 x 254	1613 x 204	550 - 880
143-1798-XX	Prone Pelvic Band Size 1	743 x 203	1283 x 153	375 - 595
143-2798-XX	Prone Pelvic Band Size 2	902 x 241	1397 x 191	460 - 730
143-3798-XX	Prone Pelvic Band Size 3	1086 x 254	1632 x 204	550 - 880

Replace XX with colour:

07 - Orange

08 - Blue

09 - Pink





Pronestander De-Rotational Pelvic Bands

Product Code	Description	Inner Fabric (mm)	Fabric Around Frame (mm)	Hip Circumference (mm)
LPS/1/B-MW	Pronestander Pelvic Band Size 1	700 x 100	750 x 100	400 – 550
LPS/2/B-MB	Pronestander Pelvic Band Size 2	900 x 100	1020 x 100	500 – 700
LPS/3/B-MB	Pronestander Pelvic Band Size 3	1360 x 100	1250 x 100	750 – 1060

N.B. These can also be used on the Leckey Freestander





Freestander De-Rotational Pelvic Bands

Product Code	Description	Inner Fabric (mm)	Fabric Around Frame (mm)	Hip Circumference (mm)
LFS/1/B-MW	Freestander Pelvic Band Size 1	640 x 100	1040 x 100	400 – 550
LFS/2/B-MB	Freestander Pelvic Band Size 2	890 x 100	1310 x 100	500 – 700
LFS/3/B-MB	Freestander Pelvic Band Size 3	1340 x 100	1530 x 100	750 – 1060

N.B. These can also be used on the Leckey Pronestander

