



Leckey BeMe



CE

Name	
Date:	
Email	
Contact Number	
Invoice Address	
Delivery Address	

Quote Order

Purchase Order	
----------------	--



	Item Code	Description	Quantity	Price
Seat Core				
				
<input type="checkbox"/>	250-1601	BeMe Size 1 Seat Core		XXX
<input type="checkbox"/>	250-2601	BeMe Size 2 Seat Core		XXX
<input type="checkbox"/>	250-3601	BeMe Size 3 Seat Core		XXX
<input type="checkbox"/>	250-4601	BeMe Size 4 Seat Core		XXX
Chassis				
				
Size 1 Indoor Chassis (Select chassis type and Push handle)				
<input type="checkbox"/>	150-1MHT-1	Multi height and tilt Indoor Manual Size 1		XXX
<input type="checkbox"/>	150-1MHT-PWR-1	Multi height and tilt Indoor Electric Size 1		XXX
<input type="checkbox"/>	150-1MHT-PH	Push Handle Size 1		XXX
Size 1 Indoor/Outdoor Chassis (Select chassis type and Push handle)				
<input type="checkbox"/>	150-1MHT-1B	Multi height and tilt Indoor-Outdoor Manual Size 1		XXX
<input type="checkbox"/>	150-1MHT-PWR-1B	Multi height and tilt Indoor-Outdoor Electric Size 1		XXX
<input type="checkbox"/>	150-1MHT-PH-1B	Push Handle Size 1 Indoor-Outdoor		XXX
Size 2 Indoor Chassis (Select chassis type and Push handle)				
<input type="checkbox"/>	150-2MHT-1	Multi height and tilt Indoor Manual Size 2		XXX
<input type="checkbox"/>	150-2MHT-PWR-1	Multi height and tilt Indoor Electric Size 2		XXX
<input type="checkbox"/>	150-2MHT-PH	Push Handle Size 2		XXX
Size 2 Indoor/Outdoor Chassis (Select chassis type and Push handle)				
<input type="checkbox"/>	150-2MHT-1B	Multi height and tilt Indoor-Outdoor Manual Size 2		XXX
<input type="checkbox"/>	150-2MHT-PWR-1B	Multi height and tilt Indoor-Outdoor Electric Size 2		XXX
<input type="checkbox"/>	150-2MHT-PH-1B	Push Handle Size 2 Indoor-Outdoor		XXX
Size 3-4 Indoor Chassis (Select chassis type and Push handle)				
<input type="checkbox"/>	150-3MHT-1	Multi height and tilt Indoor Manual Size 3		XXX
<input type="checkbox"/>	150-3MHT-PWR-1	Multi height and tilt Indoor Electric Size 3		XXX
<input type="checkbox"/>	150-3MHT-PH	Push Handle Size 3		XXX
Size 3-4 Indoor/Outdoor Chassis (Select chassis type and Push handle)				
<input type="checkbox"/>	150-3MHT-1B	Multi height and tilt Indoor-Outdoor Manual Size 3		XXX
<input type="checkbox"/>	150-3MHT-PWR-1B	Multi height and tilt Indoor-Outdoor Electric Size 3		XXX
<input type="checkbox"/>	150-3MHT-PH-1B	Push Handle Size 3 Indoor-Outdoor		XXX


Power leads (Select correct lead if ordering any Electric chassis)				
<input type="checkbox"/>	250-MHT-EUR	Adapter for EU plug		XXX
<input type="checkbox"/>	250-MHT-USA	Adapter for USA plug		XXX
<input type="checkbox"/>	250-MHT-AUS	Adapter for AUS plug		XXX
<input type="checkbox"/>	250-MHT-UK	Adapter for UK plug		XXX
Interfaces				
				
<input type="checkbox"/>	150-1760	Interface for Size 1 BeMe on to Zippie Iris, Quickie Iris, Quickie SR45		XXX
<input type="checkbox"/>	150-1771	Interface for Size 1 BeMe on to Q300 M Base -std		XXX
<input type="checkbox"/>	150-1772	Interface for Size 1 BeMe on to Q300 M Base - lift & tilt		XXX
<input type="checkbox"/>	150-2760	Interface for Size 2 BeMe on to Zippie Iris, Quickie Iris, Quickie SR45		XXX
<input type="checkbox"/>	150-2771	Interface for Size 2 BeMe on to Q300 M Base - std		XXX
<input type="checkbox"/>	150-2772	Interface for Size 2 BeMe on to Q300 M Base- lift & tilt		XXX
<input type="checkbox"/>	150-3760	Interface for Size 3-4 BeMe on to Zippie Iris, Quickie Iris or Quickie SR45		XXX
<input type="checkbox"/>	150-3771	Interface for Size 3-4 BeMe on to Q300 M Base-std		XXX
<input type="checkbox"/>	150-3772	Interface for Size 3-4 BeMe on to Q300 M Base - lift & tilt		XXX
<input type="checkbox"/>	150-773	User control mount bracket		XXX
Seat Hardware				
				
<input type="checkbox"/>	250-1637	Size 1 Basic Seat		XXX
<input type="checkbox"/>	250-1638	Size 1 Contoured Seat		XXX
<input type="checkbox"/>	250-2637	Size 2 Basic Seat		XXX
<input type="checkbox"/>	250-2639	Size 2 Split Seat		XXX
<input type="checkbox"/>	250-2638	Size 2 Contoured Seat		XXX
<input type="checkbox"/>	250-3637	Size 3 Basic Seat		XXX
<input type="checkbox"/>	250-3639	Size 3 Split Seat		XXX
<input type="checkbox"/>	250-3638	Size 3 Contoured Seat		XXX
<input type="checkbox"/>	250-4637	Size 4 Basic Seat		XXX
<input type="checkbox"/>	250-4639	Size 4 Split Seat		XXX
<input type="checkbox"/>	250-4638	Size 4 Contoured Seat		XXX
Pelvic Belt (add Bodypoint)				
				
<input type="checkbox"/>	150-1844	4-point belt Size 1		XXX
<input type="checkbox"/>	150-2844	4-point belt Size 2		XXX
<input type="checkbox"/>	150-3844	4-point belt Size 3		XXX
<input type="checkbox"/>	150-4844	4-point belt Size 4		XXX
<input type="checkbox"/>	117-476	Pelvic Cradle size 2		XXX
<input type="checkbox"/>	117-477	Pelvic Cradle size 3		XXX
<input type="checkbox"/>	117-478	Pelvic Cradle size 4		XXX
<input type="checkbox"/>	117-479	Pelvic Cradle size 5		XXX
<input type="checkbox"/>	117-480	Pelvic Cradle size 6		XXX
<input type="checkbox"/>	150-1868	Small Pelvic Cradle Attachment Strap		XXX
<input type="checkbox"/>	150-2868	Medium Pelvic Cradle Attachment Strap		XXX
<input type="checkbox"/>	150-3868	Large Pelvic Cradle Attachment Strap		XXX

Femoral Guide



<input type="checkbox"/>	150-1615	Single Femoral Gable - Small		XXX
<input type="checkbox"/>	150-2615	Single Femoral Gable - Medium		XXX
<input type="checkbox"/>	150-3615	Single Femoral Gable - Large		XXX
<input type="checkbox"/>	150-4615	Single Femoral Gable - Extra Large		XXX
<input type="checkbox"/>	150-1862	Padded Cover for Single Femoral Gable - Small		XXX
<input type="checkbox"/>	150-2862	Padded Cover for Single Femoral Gable - Medium		XXX
<input type="checkbox"/>	150-3862	Padded Cover for Single Femoral Gable - Large		XXX
<input type="checkbox"/>	150-4862	Padded Cover for Single Femoral Gable - Extra Large		XXX
<input type="checkbox"/>	150-1862-10	Padded Cover for Single Femoral Gable - Small (incon)		XXX
<input type="checkbox"/>	150-2862-10	Padded Cover for Single Femoral Gable - Medium (incon)		XXX
<input type="checkbox"/>	150-3862-10	Padded Cover for Single Femoral Gable - Large (incon)		XXX
<input type="checkbox"/>	150-4862-10	Padded Cover for Single Femoral Gable - Extra Large (incon)		XXX

Pommel








<input type="checkbox"/>	150-1624	Pommel Small		XXX
<input type="checkbox"/>	150-1864	Padded Cover - Pommel Small		XXX
<input type="checkbox"/>	150-1864-10	Padded Cover - Pommel Small (incon)		XXX
<input type="checkbox"/>	150-2624	Pommel Medium		XXX
<input type="checkbox"/>	150-2864	Padded Cover - Pommel Medium		XXX
<input type="checkbox"/>	150-2864-10	Padded Cover - Pommel Medium (incon)		XXX
<input type="checkbox"/>	150-3624	Pommel Large		XXX
<input type="checkbox"/>	150-4624	Pommel X Large		XXX
<input type="checkbox"/>	150-3864	Padded Cover - Pommel Large		XXX
<input type="checkbox"/>	150-3864-10	Padded Cover - Pommel Large (incon)		XXX


Backrest



<input type="checkbox"/>	250-1625	Size 1 Basic Backrest		XXX
<input type="checkbox"/>	250-1630	Size 1 Moderate Backrest		XXX
<input type="checkbox"/>	250-2625	Size 2 Basic Backrest		XXX
<input type="checkbox"/>	250-2630	Size 2 Moderate Backrest		XXX
<input type="checkbox"/>	250-2631	Size 2 Complex 3D Backrest		XXX
<input type="checkbox"/>	250-3625	Size 3 Basic Backrest		XXX
<input type="checkbox"/>	250-3630	Size 3 Moderate Backrest		XXX
<input type="checkbox"/>	250-3631	Size 3 Complex 3D Backrest		XXX
<input type="checkbox"/>	250-4625	Size 4 Basic Backrest		XXX
<input type="checkbox"/>	250-4630	Size 4 Moderate Backrest		XXX
<input type="checkbox"/>	250-4631	Size 4 Complex 3D Backrest		XXX
<input type="checkbox"/>	150-1792	15 degree angle spacer for Moderate shoulder section		XXX
<input type="checkbox"/>	150-605	Dynamics Upgrade Kit		XXX

Standard Hip Guides				
				
<input type="checkbox"/>	150-1684	Small Hip Guide (Pair)		XXX
<input type="checkbox"/>	150-2684	Medium Hip Guide (Pair)		XXX
<input type="checkbox"/>	150-3684	Large Hip Guide (Pair)		XXX
<input type="checkbox"/>	150-1684-10	Small Hip Guide (Pair) (incon)		XXX
<input type="checkbox"/>	150-2684-10	Medium Hip Guide (Pair) (incon)		XXX
<input type="checkbox"/>	150-3684-10	Large Hip Guide (Pair) (incon)		XXX
Slim Hip Guide				
				
<input type="checkbox"/>	150-1678	Small Slim Hip Guide (Pair)		XXX
<input type="checkbox"/>	150-2678	Medium Slim Hip Guide (Pair)		XXX
<input type="checkbox"/>	150-3678	Large Slim Hip Guide (Pair)		XXX
<input type="checkbox"/>	150-1678-10	Small Slim Hip Guide (Pair) (incon)		XXX
<input type="checkbox"/>	150-2678-10	Medium Slim Hip Guide (Pair) (incon)		XXX
<input type="checkbox"/>	150-3678-10	Large Slim Hip Guide (Pair) (incon)		XXX
<input type="checkbox"/>	250-1682	Inset small Hip Guide (Pair)		XXX
Hip Guides - Additional Longer Pads (must order hip guides)				
				
<input type="checkbox"/>	150-1653	Hinged Lateral Pad Small (Pair)		XXX
<input type="checkbox"/>	150-3653	Hinged Lateral Pad Large (Pair)		XXX
<input type="checkbox"/>	150-1653-10	Hinged Lateral Pad Small (Pair) (incon)		XXX
<input type="checkbox"/>	150-3653-10	Hinged Lateral Pad Large (Pair) (incon)		XXX
Chest Harness				
				
<input type="checkbox"/>	155-N001-L600	Small Dynamic Butterfly Harness		XXX
<input type="checkbox"/>	155-N001-L601	Medium Dynamic Butterfly Harness		XXX
<input type="checkbox"/>	155-N001-L602	Large Dynamic Butterfly Harness		XXX
<input type="checkbox"/>	155-N001-L603	Extra Large Dynamic Butterfly Harness		XXX
Lateral Bracket				
				
<input type="checkbox"/>	250-1651	Rigid Lateral Bracket for Size 1/2 Backrests(Pair)		XXX
<input type="checkbox"/>	250-1655	FlipAway Lateral Bracket for Size 1/2 Backrests (Pair)		XXX
<input type="checkbox"/>	250-3651	Rigid Lateral Bracket for Size 3/4 Backrests (Pair)		XXX
<input type="checkbox"/>	250-3655	FlipAway Lateral Bracket for Size 3/4 Backrests (Pair)		XXX
<input type="checkbox"/>	250-1661	Inset Lateral Bracket (Pair)		XXX
<input type="checkbox"/>	250-1665	Inset Flipaway Lateral Bracket (Pair)		XXX

Chest Lateral



<input type="checkbox"/>	150-1675	Size 1 Slimline Lateral Pads (Pair)		XXX
<input type="checkbox"/>	150-2675	Size 2 Slimline Lateral Pads (Pair)		XXX
<input type="checkbox"/>	150-3675	Size 3 Slimline Lateral Pads (Pair)		XXX
<input type="checkbox"/>	150-1260	Small Wraparound Hugging Lateral Pads (Pair)		XXX
<input type="checkbox"/>	150-2260	Medium Wraparound Hugging Lateral Pads (Pair)		XXX
<input type="checkbox"/>	150-1653	Hinged Lateral Pad Small (Pair)		XXX
<input type="checkbox"/>	150-3653	Hinged Lateral Pad Large (Pair)		XXX
<input type="checkbox"/>	150-673	Lateral Protraction Pads (Pair)		XXX
<input type="checkbox"/>	150-1675-10	Size 1 Slimline Lateral Pads (Pair) (incon)		XXX
<input type="checkbox"/>	150-2675-10	Size 2 Slimline Lateral Pads (Pair) (incon)		XXX
<input type="checkbox"/>	150-3675-10	Size 3 Slimline Lateral Pads (Pair) (incon)		XXX
<input type="checkbox"/>	150-1260-10	Small Wraparound Hugging Lateral Pads (Pair) (incon)		XXX
<input type="checkbox"/>	150-2260-10	Medium Wraparound Hugging Lateral Pads (Pair) (incon)		XXX
<input type="checkbox"/>	150-1653-10	Hinged Lateral Pad Small (Pair) (incon)		XXX
<input type="checkbox"/>	150-3653-10	Hinged Lateral Pad Large (Pair) (incon)		XXX
<input type="checkbox"/>	150-673-10	Lateral Protraction Pads (Pair) (incon)		XXX

Armrest



<input type="checkbox"/>	150-1642	Small Armrests (Pair)		XXX
<input type="checkbox"/>	150-2642	Medium Armrests (Pair)		XXX
<input type="checkbox"/>	150-3642	Large Armrests (Pair)		XXX

Additional Armrest Pad (single)



<input type="checkbox"/>	150-15204	Small Armrest Pad Only (Single)		XXX
<input type="checkbox"/>	150-25204	Medium/Large Armrest Pad Only (Single)		XXX

Tray (match Seat core size)



<input type="checkbox"/>	150-1720	Size 1 Tray		XXX
<input type="checkbox"/>	150-1888	Size 1 Padded Tray Cover		XXX
<input type="checkbox"/>	150-1671	Tray Elbow Blockers for Size 1		XXX
<input type="checkbox"/>	150-2720	Size 2 Tray		XXX
<input type="checkbox"/>	150-2888	Size 2 Padded Tray Cover		XXX
<input type="checkbox"/>	117-641	Tray Elbow Blockers for Size 2		XXX
<input type="checkbox"/>	150-3720	Size 3 Tray		XXX
<input type="checkbox"/>	150-3888	Size 3 Padded Tray Cover		XXX
<input type="checkbox"/>	150-3671	Tray Elbow Blockers for Size 3		XXX
<input type="checkbox"/>	150-4720	Size 4 Tray	(Use Size 3 cover & Elbow Blockers if required)	XXX

Leg Support				
				
<input type="checkbox"/>	250-1689	Split Footplate Size 1(Left and Right)		XXX
<input type="checkbox"/>	250-2689	Split Footplate Size 2 (Left and Right)		XXX
<input type="checkbox"/>	250-3689	Split Footplate Size 3 (Left and Right)		XXX
<input type="checkbox"/>	250-4689	Split Footplate Size 4 (Left and Right)		XXX
<input type="checkbox"/>	117-787	1Inch Sandal Raisers (pair) for split footplate		XXX
<input type="checkbox"/>	250-1698	One Piece Footplate - Size 1		XXX
<input type="checkbox"/>	250-2698	One Piece Footplate - Size 2		XXX
<input type="checkbox"/>	250-3698	One Piece Footplate - Size 3		XXX
<input type="checkbox"/>	250-4698	One Piece Footplate - Size 4		XXX
<input type="checkbox"/>	150-1872	Cushion for One-Piece Footplate - Size 1		XXX
<input type="checkbox"/>	150-2872	Cushion for One-Piece Footplate - Size 2		XXX
<input type="checkbox"/>	150-3872	Cushion for One-Piece Footplate - Size 3/4		XXX
Foot Support				
				
<input type="checkbox"/>	AH01	Ankle Huggers		XXX
<input type="checkbox"/>	150-1852	Extra Small Booties		XXX
<input type="checkbox"/>	150-2852	Small Booties		XXX
<input type="checkbox"/>	152-1600	Sandals Small		XXX
<input type="checkbox"/>	152-2600	Sandals Medium		XXX
<input type="checkbox"/>	152-3600	Sandals Large		XXX
<input type="checkbox"/>	173-2305	Small Ankle Wraps (Only with Sandals)		XXX
<input type="checkbox"/>	173-3305	Large Ankle Wraps (Only with Sandals)		XXX
Headrests				
				
<input type="checkbox"/>	150-675	BeMe Headrest		XXX
<input type="checkbox"/>	150-780	Universal Headrest Bracket (Whitmyer PRO or Leckey Headrests)		XXX
Leckey Headrests				
<input type="checkbox"/>	137-699-05	Leckey Contoured Headrest		XXX
<input type="checkbox"/>	117-699-05	Leckey Contoured Occipital Headrest		XXX
<input type="checkbox"/>	117-851-05	Leckey Flat Headrest		XXX
<input type="checkbox"/>	120-893-05	Leckey Flat Headrest Laterals		XXX
Miscellaneous				
<input type="checkbox"/>	150-1745	Communication mount		XXX
<input type="checkbox"/>	150-1645	Sun/Rain Canopy for Small Backrest		XXX
<input type="checkbox"/>	150-2645	Sun/Rain Canopy for Medium Backrest		XXX
<input type="checkbox"/>	150-1790	Spacer Pair to decrease minimum seat depth on BeMe Small - 25mm		XXX
<input type="checkbox"/>	150-1791	Spacer Pair to decrease minimum seat depth on BeMe Small - 50mm		XXX
<input type="checkbox"/>	150-846	Extra Padded Sleeve for Footplate Stem or Lateral Brackets		XXX

Upholstery				
		Blue - Replace XX with 02		
		Purple - Replace XX with 04		
		Red - Replace XX with 09		
		Grey - Replace XX with 08		
		Incon - Replace XX with 10		
<input type="checkbox"/>	150-855-XX	BeMe Headrest Cushion and Cover		XXX
<input type="checkbox"/>	250-1840-XX	Size 1 Basic Seat Cushion and Cover		XXX
<input type="checkbox"/>	250-1860-XX	Size 1 Contoured Seat Cushion and Cover		XXX
<input type="checkbox"/>	250-2840-XX	Size 2 Basic Seat Cushion and Cover		XXX
<input type="checkbox"/>	250-2860-XX	Size 2 Contoured Seat Cushion and Cover		XXX
<input type="checkbox"/>	250-2850-XX	Size 2 Split Seat Cushion and Cover		XXX
<input type="checkbox"/>	250-3840-XX	Size 3 Basic Seat Cushion and Cover		XXX
<input type="checkbox"/>	250-3860-XX	Size 3 Contoured Seat Cushion and Cover		XXX
<input type="checkbox"/>	250-3850-XX	Size 3 Split Seat Cushion and Cover		XXX
<input type="checkbox"/>	250-4840-XX	Size 4 Basic Seat Cushion and Cover		XXX
<input type="checkbox"/>	250-4860-XX	Size 4 Contoured Seat Cushion and Cover		XXX
<input type="checkbox"/>	250-4850-XX	Size 4 Split Seat Cushion and Cover		XXX
<input type="checkbox"/>	250-1870-XX	Size 1 Basic Backrest Cover		XXX
<input type="checkbox"/>	250-1800-XX	Size 1 Moderate Backrest Cushions and Covers		XXX
<input type="checkbox"/>	250-2870-XX	Size 2 Basic Backrest Cover		XXX
<input type="checkbox"/>	250-2800-XX	Size 2 Moderate Backrest Cushions and Covers		XXX
<input type="checkbox"/>	250-4800-XX	Size 2 Complex Backrest Cushions and Covers		XXX
<input type="checkbox"/>	250-3870-XX	Size 3 Basic Backrest Cover		XXX
<input type="checkbox"/>	250-3800-XX	Size 3 Moderate Backrest Cushions and Covers		XXX
<input type="checkbox"/>	250-5800-XX	Size 3 Complex Backrest Cushions and Covers		XXX
<input type="checkbox"/>	250-4870-XX	Size 4 Basic Backrest Cover		XXX
<input type="checkbox"/>	250-6800-XX	Size 4 Moderate Backrest Cushions and Covers		XXX
<input type="checkbox"/>	250-7800-XX	Size 4 Complex Backrest Cushions and Covers		XXX

Technical Information

	Size 1	Size 2	Size 3	Size 4
Seat Depth (mm) *	120-270	225-350	320-430	400-540
Seat Width Maximum (No Hip Guides) (mm)	280	325	410	475
Femoral Width on Split Seat (mm)		90-125	110-160	135-170
Maximum Weight Limit (kg) **	35	50	70	80
	Size 1	Size 2	Size 3	Size 4
Back Height (mm)	300-400	350-470	400-530	475-650
Standard Hip Guide Width (mm)	145-215	175-260	245-350	290-405
Slim Hip Guide Width (mm)	170-240	200-285	270-375	315-430
Lateral Width (mm)	130-230	200-300	225-350	350-450
Inset Lateral Width (mm)	80-180	150-250	175-300	300-400
Backrest Angle (degrees)			+10 to -25	
Maximum Weight Limit (kg) **	35	50	70	80